



## MODEL MM1P

THIS 100 GALLON MIXING TANK, EMPLOYS JET FLOW INJECTION FOR MIXING WITH A SIMPLE ON/OFF SWITCH. WE USE A MIXING PUMP WITH THE POTENTIAL OF 53 GALLONS PER MINUTE UNRESTRICTED FLOW TO OUR MIXING JET TO INSURE THE CONCENTRATE IS MIXED INSTANTLY WITH A TURBULANCE MULTIPLIER FORCE IN EXCESS OF A ONE HUNDRED GALLON PER MINUTE OUT FLOW.

## MODEL MM5P

THIS 50 GALLON MIXING TANK, EMPLOYS JET FLOW INJECTION FOR MIXING WITH A SIMPLE ON/OFF SWITCH. WE USE A MIXING PUMP WITH THE POTENTIAL OF 53 GALLONS PER MINUTE UNRESTRICTED FLOW TO OUR MIXING JET TO INSURE THE CONCENTRATE IS MIXED INSTANTLY WITH A TURBULANCE MULTIPLIER FORCE IN EXCESS OF A ONE HUNDRED GALLON PER MINUTE OUT FLOW.

## AUTO FILL

OPTIONAL AUTO FILL CONTROLLER IS AVAILABLE FOR BOTH MODELS TO ALLOW FOR UNATTENDED FILLING OF THE TANK.

MEDICAL SOLUTIONS INTERNATIONAL, INC.

11522 WEST 90TH STREET, SHAWNEE MISSION, KS 66214  
PH: 913-438-9700 FAX: 913-438-9701

UNIT: MM1P	TITLE: 100 GAL. BICARB MIXER		
MATERIAL: N/A	SCALE: NONE	APPROVED BY: JACK R. DILLON	DRN: JLS
FINISH: N/A	DATE: 05/19/04	NEXT ASSEMBLY	REV: 1
		DOCUMENT #	SHT 1
			OF 1

Medical Solutions International, Inc. (MSI) PROPRIETARY  
This document and all information contained herein are the property of MSI, Inc., are provided in confidence, and may not be used, duplicated or disclosed to a third party, in whole or in part, without prior written permission of MSI, Inc. Copyright 2003 Medical Solutions Int., Inc. All rights reserved.

TOL UNLESS NOTED:  
.XX ±  
.XXX ±  
ANGLES ±