


COILED CONCENTRATE LINES



Use for connections between your solutions delivery equipment and your dialysis machines. We recommend our Coily's as the perfect solution for interfacing with both our standard and recessed Patient Wall Stations.

Each Coily is translucent for ease of fluid verification, and may be color-coordinated to your system.

Polyurethane has outstanding memory and flexibility which makes it ideal material for processing into self-storing Coilys. The coil allows the line to retract off the floor as equipment is moved, keeping the area neat. These coils are extraordinarily tough, resistant to abrasion, overstretching and kniking. They are very flexible and come with two straight tails, allowing for attachment at a comfortable distance from your equipment.

Standard Coilys are 10' long when fully extended, 1/4" OD and 1/8" ID. They are available in the following colors:

PART NO.	COLOR	PRICE
C80550	RED (Translucent)	\$11.95
C80553	BLUE (Translucent)	\$11.95
C80558	CLEAR	\$11.95
C80558-B	CLEAR - for Baxter 1550	\$11.95

TO ORDER

CALL 1-800-326-5275

MEDICAL SOLUTIONS INTERNATIONAL, INC.
 11522 WEST 90TH STREET, SHAWNEE MISSION, KS 66214
 PH: 913-438-9700 FAX: 913-438-9701

PATENT APPLIED FOR

Medical Solutions International, Inc. (MSI) PROPRIETARY
 This document and all information contained herein are the property of MSI, Inc., are provided in confidence, and may not be used, duplicated or disclosed to a third party, in whole or in part, without prior written permission of MSI, Inc. Copyright 2003 Medical Solutions Int., Inc. All rights reserved.

TOL UNLESS NOTED:
 .XX ±
 .XXX ±
 ANGLES ±

UNIT: —	TITLE: Coily		
MATERIAL: N/A	SCALE: NONE	APPROVED BY: JACK R. DILLON	DRN: JLS
FINISH: N/A	DATE: 05/18/04	NEXT ASSEMBLY	REV: 1
		DOCUMENT #	SHT OF
		N/A	— OF