



Product Catalog

2022





Medical Solutions International has been privileged to serve the dialysis community since 1989. We designed and manufactured the first commercial Bicarbonate Mixing and Distribution System, the first commercial Acidified Distribution systems and the first dialysis specific commercial Patient Wall Station.

We have built our business on the faith and trust of our dialysis customers. Our innovative equipment designs were improved over the years by listening to our customers needs and performance expectations.

We will continue to strive to be innovative and reactive to our customers. Please feel free to contact us if you have an idea that you feel could be developed to better serve the dialysis community.



Medical Solutions

I N T E R N A T I O N A L



TABLE OF CONTENTS

Call: 1-800-326-5275
Available Monday thru
Friday 8:00AM-4:00PM CST

Email: sales@medicalsolutionskc.com
Customer support available
Monday thru Friday 8:00AM-4:00pm CST

Catalog also available on-line
www.medicalsolutionskc.com

SDS-CB200 DUAL TANK BICARB SYSTEM	1
NEW PRODUCTS -RACK, BLEACH CAP, ETC	2
CITRIC JUG TEE, JUG CAP & JUG TAP	3
ACID CABINET & CHLORINE TIMER	4
MR. ACID	5
CRANK HOLDER/SAMPLE CENTER	6
OZONE ACCESSORIES	7
WALL TECH STATION	8-12
WALL TECH STATION ADAPTERS	13-16
CLEAR SCREEN PROTECTOR	17
BICARB SYSTEM ACCESSORIES	18-20
FUSES / SWITCHES / SENSORS	21
ELECTRICAL	22
PRESSURIZED PRODUCTS	23
CAMLOCKS / BANJO BULKHEAD	24
FILTERS	25
PUMPS	26-27
HOSE / CLAMPS / CLICS	28
TUBING	29
PVC SCHEDULE 80 FITTINGS	30-41
NIPPLES	42-44
BALL VALVES	45
CHECK VALVES / GATE VALVES	46
BARBED FITTINGS	47
QUICK DISCONNECT FITTINGS	48-59
PARKER FITTINGS	60-62
MUR-LOK	63-66
GARDEN HOSE FITTINGS	67
TANKS	68
ORDER FORM	69

Medical Solutions

INTERNATIONAL

Bicarbonate Mixing & Distribution System




Dual Tank Bicarbonate Mixing and Pressurized Distribution System

These systems have 25 to 100 gallon mixing capacity, with electronic fill and dispense options, variable batch size, internally balanced agitation, sanitization/disinfection mode, a .2 micron air filter, and timed cycle processing, they employ two 100 gallon tanks, one for mixing and one for distribution. With this system, you may effectively reconstitute dialysate powders in a closed system with a managed fill to accomplish predictable/reproducible batch quality.

The systems are designed for low volume, low pressure dynamic distribution of mixed Bicarb concentrate to patient stations. The system lends them self to simple, complete sanitization and effective rapid rinse. The system is

SDS-CB200 DUAL TANK BICARBONATE MIXING & DISTRIBUTION 1



CALL FOR BICARBONATE SYSTEMS QUOTATIONS

ELECTRICAL CONTROL BOXES


PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
M-RA02	Bicarb Remote Alarm with Reset Mute	641.25



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

Wand Drying Rack




CBM-DR-1
Wand Drying Rack



medicalsolutionskc.com 800-326-5275

DRAIN HOSE CONTROL FOR ACUTE TREATMENT



medicalsolutionskc.com 800-326-5275

Bleach Cap

Cap Only
(Jar not included)



medicalsolutionskc.com 800-326-5275

NEW PRODUCTS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CBM-DR-1	Wand Drying Rack	29.36
DHH-1	Butterfly Hose Holder	12.90
K-FRE-BC	Bleach Cap-All Plastic	26.60
K-FRE-BOTTLE	500 ml Bottle (cap not included)	13.27

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701


JUG CAPS, TAPS & TEES

3

Red Jug Cap



M-JC-R
FOR FRESenius
ACIDIFIED CONNECTORS




M-JC-R-BB
FOR B.BRAUN
ACIDIFIED CONNECTORS

Blue Jug Cap



M-JC-B
FOR FRESenius
BICARBONATE CONNECTORS



M-JC-B-BB
FOR B.BRAUN
BICARBONATE CONNECTORS

"Clean, Neat, No Mess".



medicalsolutionskc.com 800-326-5275



medicalsolutionskc.com 800-326-5275

Red Jug Tap



M-JT-R

Blue Jug Tap



M-JT-B

"Clean, Neat, No Mess".

JUG CAP, TAP & TEES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CAT-3	Citric Acid Jug Tee	49.20
CAT-2	Vinegar Jug Tee	49.20
M-JC-R	Red Acid Jug Cap (Fresenius)	26.60
M-JC-B	Blue Bicarb Jug Cap (Fresenius)	26.60
M-JC-R-BB	Red Acid Cap (BBraun)	26.60
M-JC-B-BB	Blue Bicarb Cap (BBraun)	26.60
M-JT-R	Red Acid Jug Tap (Fresenius)	22.10
M-JT-B	Blue Bicarb Jug Tap (Fresenius)	22.10

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

MEDICAL SOLUTIONS INTERNATIONAL




ACIDIFIED CONCENTRATE DISTRIBUTION SYSTEM

We use all corrosion resistant materials as well as FDA registered process piping, tanks, controls, etc. The system is designed to transfer acidified concentrate from bulk storage to the patient stations by controlled low pressure pump. The system consists of transfer pumps, valves, fittings, and a complete electronic control. Includes pressurized jug fill.


"It won't let you forget".

Chlorine timer



M-CT02

Remote Alarm



M-RA03



medicalsolutionskc.com 800-326-5275

ACIDIFIED DISTRIBUTION AND CHLORINE TIMER/REMOTE

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
APDS-CAB3	Acidified Distribution Cabinet	2954.00
M-CT02	Chlorine Timer	586.95
M-RA03	Remote Alarm for Chlorine Timer	129.95

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701



MR. ACID

If They Want You To Fill Jugs The Old Way...
JUST SAY NO!

*This Is Just Too Easy!
Automatic Jug Filling!
Clean, Spill Free, Easy
Shuts Itself Off When Full...*

Just place the jug, push the button, walk away...
Mr. Acid! Will fill the jug and automatically shut off
when full. Mr. Acid! Works on any dialysis machine
jug and works with any acidified concentrate or with
bicarb.

Mr. Acid! Is clean, spill free and easy. No
more slick dangerous floors. No more
spilling and sliding. No more ruined shoes
and floors.

This system is made to fit right on the acid drum
or can be ordered for carts to fill from bulk
storage tanks.

This Patent Applied For



ADU-55



medicalsolutionskc.com 800-326-5275

MR. ACID		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
M-ADU-55	Mr. Acid	1318.45
M-A-DW01	Mr. Acid Drum Wand	21.53



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com



800.326.5275
 913.438.9700
 FAX 913.438.9701

SAMPLE— —CENTER

- EASY
- SIMPLE
- CONVENIENT
- CLEAN

STORE YOUR TESTING EQUIPMENT
CONVENIENTLY AWAY...

Medical Solutions
INTERNATIONAL
www.medicalsolutionskc.com
800-326-5275



Medical Solutions
INTERNATIONAL

medicalsolutionskc.com 800-326-5275

SAMPLE STATION

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
M-FCH-001	Machine Crank Holder	11.94
SS-101	Pre-Treatment Sample Station w/out Gauges	533.25
SS-101-PG	Pre-Treatment Sample Station w/ Gauges	845.25

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX **913.438.9701**



Z-C30ZX



O3-DRF-8



BEFORE AFTER



O3-T

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
O3-DRF-H1	Ozone Dessicant Filter Housing	39.95
O3-DRF-8	Ozone Air Drying Dessicant Refill	11.91
Z-C-30ZX	Ozone Ambient Air Monitor	695.75
Z-CX01	Exchange Cell for Ozone Unit	269.98
O3-T	Test Strips for Ozone (1 Bottle)	29.95
CV-14000	1/4" Check Valve	18.83

MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com


800.326.5275
 913.438.9700
 FAX 913.438.9701

Drain Isolation Wall Station



Slide
Door Open

NEW



CALL FOR PRICING ON NEW WALLSTATION,
MULTIPLE CONFIGURATIONS TO MEET YOUR NEEDS!

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

Medical Solutions International




Wall Station for Dialysis

Our recessed dialysis wall stations include a solutions distribution center and drain. Each is 18" wide by 18" high, and made of high quality molded polystyrene. The durable material can be fitted with a variety of valves, fittings, pressure gauges and quick connects so as to accommodate the installation of compatible dispensing with any dialysis delivery system on the market today. The easily cleaned, chemical resistant wall station also receives the drainage from the dialysis machine. The drip-proof canister valve block keeps the unit free of spills and stains. All valving is color coded and signage designated.

- A. P6FC8
- B. ML-T0620666
- C. ML-G0620625
- D. W-M11X4 (4 SLOT MANIFOLD)
- E. W-FRE-850 (RED, BLUE, YELLOW, WHITE)
- F. QD-HFCD24812
- G. QD-HFCD101235GHT
- H. HM001(BLUE, RED, GREEN, ORANGE, YELLOW)
- I. W-CNC-AGD (AIR GAP DRAIN ASSEMBLY)

NOTE:
TUBING (NOT SHOWN) SHALL BE INSTALLED
AS TO THE TYPICAL 4-PLACEMENT.



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701


WALL TECH STATIONS-BEIGE

10

BEIGE WALL STATIONS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-RPWS-B000	Beige Wallstation Blank	214.65
W-RPWS-B012	Beige Wallstation with 1 Bicarb and 2 Acid Fittings	389.14
W-RPWS-B100	Beige Wallstation with 1 Water Fitting	282.08
W-RPWS-B101	Beige Wallstation with 1 Water and 1 Acid Fitting	336.96
W-RPWS-B102	Beige Wallstation with 1 Water and 2 Acid Fittings	389.14
W-RPWS-B103	Beige Wallstation with 1 Water and 3 Acid Fittings	443.88
W-RPWS-B104	Beige Wallstation with 1 Water and 4 Acid Fittings	497.81
W-RPWS-B111	Beige Wallstation with 1 Water, 1 Bicarb and 1 Acid Fitting	389.95
W-RPWS-B112	Beige Wallstation with 1 Water, 1 Bicarb and 2 Acid Fittings	443.88
W-RPWS-B113	Beige Wallstation with 1 Water, 1 Bicarb and 3 Acid Fittings	497.81
W-RPWS-B114	Beige Wallstation with 1 Water, 1 Bicarb and 4 Acid Fittings	551.75
W-RPWS-B200	Beige Wallstation with 2 Water Fittings	350.86
W-RPWS-B202	Beige Wallstation with 2 Water and 2 Acid Fittings	457.38
W-RPWS-B212	Beige Wallstation with 2 Water, 1 Bicarb and 2 Acid Fittings	478.58
W-RPWS-B222	Beige Wallstation with 2 Water, 2 Bicarb and 2 Acid Fittings	565.25
W-RPWS-B300	Beige Wallstation with 3 Water Fittings	416.95
W-RPWS-B400	Beige Wallstation with 4 Water Fittings	443.88

REFER TO PAGE 9



NOTE:
TURNING (NOT SHOWN) PLATE FROM
AS TO OIL TYPICAL PLACES

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701


WALL TECH STATIONS-GRAY

11

GRAY WALL STATIONS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-RPWS-G000	Gray Wallstation Blank	214.65
W-RPWS-G012	Gray Wallstation with 1 Bicarb and 2 Acid Fittings	389.14
W-RPWS-G100	Gray Wallstation with 1 Water Fitting	282.08
W-RPWS-G101	Gray Wallstation with 1 Water and 1 Acid Fitting	336.96
W-RPWS-G102	Gray Wallstation with 1 Water and 2 Acid Fittings	389.95
W-RPWS-G103	Gray Wallstation with 1 Water and 3 Acid Fittings	443.88
W-RPWS-G104	Gray Wallstation with 1 Water and 4 Acid Fittings	497.81
W-RPWS-G111	Gray Wallstation with 1 Water, 1 Bicarb and 1 Acid Fitting	389.95
W-RPWS-G112	Gray Wallstation with 1 Water, 1 Bicarb and 2 Acid Fittings	443.88
W-RPWS-G113	Gray Wallstation with 1 Water, 1 Bicarb and 3 Acid Fittings	497.81
W-RPWS-G114	Gray Wallstation with 1 Water, 1 Bicarb and 4 Acid Fittings	551.75
W-RPWS-G200	Gray Wallstation with 2 Water Fittings	350.86
W-RPWS-G202	Gray Wallstation with 2 Water and 2 Acid Fittings	457.38
W-RPWS-G212	Gray Wallstation with 2 Water, 1 Bicarb and 2 Acid Fittings	478.58
W-RPWS-G222	Gray Wallstation with 2 Water, 2 Bicarb and 2 Acid Fittings	565.24
W-RPWS-G300	Gray Wallstation with 3 Water Fittings	416.95
W-RPWS-G400	Gray Wallstation with 4 Water Fittings	443.88

REFER TO PAGE 9



NOTE:
TIGHTEN FIRST ON KNOBS, THEN
ON TO ALL TYPICAL PLACES.

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

WHITE WALL STATIONS


PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-RPWS-W000	White Wallstation Blank	214.65
W-RPWS-W100	White Wallstation with 1 Water Fitting	282.95
W-RPWS-W101	White Wallstation with 1 Water and 1 Acid Fitting	336.96
W-RPWS-W102	White Wallstation with 1 Water and 2 Acid Fittings	389.95
W-RPWS-W103	White Wallstation with 1 Water and 3 Acid Fittings	443.88
W-RPWS-W111	White Wallstation with 1 Water, 1 Bicarb and 1 Acid Fitting	389.95
W-RPWS-W112	White Wallstation with 1 Water, 1 Bicarb and 2 Acid Fittings	443.88
W-RPWS-W113	White Wallstation with 1 Water, 1 Bicarb and 3 Acid Fittings	497.81
W-RPWS-W114	White Wallstation with 1 Water, 1 Bicarb and 4 Acid Fittings	551.75
W-RPWS-W200	White Wallstation with 2 Water Fittings	349.52
W-RPWS-W202	White Wallstation with 2 Water and 2 Acid Fittings	457.38
W-RPWS-W212	White Wallstation with 2 Water, 1 Bicarb and 2 Acid Fittings	478.58
W-RPWS-W222	White Wallstation with 2 Water, 2 Bicarb and 2 Acid Fittings	565.24
W-RPWS-W300	Gray Wallstation with 3 Water Fittings	416.95
W-RPWS-W400	Gray Wallstation with 4 Water Fittings	443.88

SMALL WHITE WALL STATIONS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-RPWS-SW000	Small White Wallstation Blank	214.65
W-RPWS-SW00AG	Small White Wallstation with Air Gap Drain Assembly	244.60
W-RPWS-SW100	Small White Wallstation with 1 Water Fitting	282.08
W-RPWS-SW102	Small White Wallstation with 1 Water and 2 Acid Fittings	389.95
W-RPWS-SW110	Small White Wallstation with 1 Water and 1 Bicarb Fitting	336.96
W-RPWS-SW111	Small White Wallstation with 1 Water, 1 Bicarb and 1 Acid Fitting	389.95
W-RPWS-SW112	Small White Wallstation with 1 Water, 1 Bicarb and 2 Acid Fittings	443.88
W-RPWS-SW200	Small White Wallstation with 2 Water Fittings	348.97

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701



ADAPTERS		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-FRE-R-850	Red 1/2" Male NPT X Fresenius Reducer Bushing	16.13
W-FRE-B-850	Blue 1/2" Male NPT X Fresenius Reducer Bushing	16.13
W-FRE-W-850	White 1/2" Male NPT X Fresenius Reducer Bushing	16.13
W-FRE-AP-R	Red All Plastic Checkless Fresenius Reducer Bushing	22.46
W-FRE-AP-B	Blue All Plastic Checkless Fresenius Reducer Bushing	22.46
W-FRE-NC-R	Red Checkless Male NPT X Fresenius Reducer Bushing w/SS Ring	29.95
W-FRE-NC-B	Blue Checkless Male NPT X Fresenius Reducer Bushing w/SS Ring	29.95
W-PHX-RP-R	Red Phoenix Rinse Port	29.95
W-PHX-RP-B	Blue Phoenix Rinse Port	29.95
W-PHX-RP-W	White Phoenix Rinse Port	29.95
W-FRE-RS	Rubber Seal for Fresenius Rinse Port for W-FRE-850's	1.10

MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

WALL TECH STATION ADAPTERS

14

ITEMS FOR ADAPTING FRESENIUS FITTINGS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-FRE-W5-B	Fresenius Assembly Washer Blue	0.38
W-FRE-W5-R	Fresenius Assembly Washer Red	0.38
W-FRE-W5-W	Fresenius Assembly Washer White	0.38
W-FRE-W5-Y	Fresenius Assembly Washer Yellow	0.38
W-FRE-W.50	Fresenius Bushing Spacer White	0.38



WALL STATION VALVE HANDLES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
HB001	Blue Handle SMC	1.69
HR001	Red Handle SMC	1.69
HC001	Clear Handle SMC	1.69
HG001	Gray Handle SMC	1.69
HBL001	Black Handle SMC	1.69
HMB001	Blue Handle Murlock	1.69
HMR001	Red Handle Murlock	1.69
HMO001	Orange Handle Murlock	1.69
HMY001	Yellow Handle Murlock	1.69
HMGR001	Green Handle Murlock	1.69
HPY001	Yellow Parker Ballvalve Handle	3.94



HANDLES IN MURDOCK GREEN



WALL STATION NYLON DRAIN BARBS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-ATE2	Wallstation Drain 1" MNPT X 3/4" Barb Elbow	3.70
W-ATA1	Wallstation Drain 1" MNPT X 3/4" Barb Straight	3.10
436-133	Wallstation Drain 1" MNPT X 1-1/2" Socket Schedule 40	4.95

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

WALL TECH STATION ITEMS

15



MEDICAL SOLUTIONS INTERNATIONAL 800-326-5275

35X BICARB	H ₂ O WATER	35X ACIDIFIED	35X ACIDIFIED	35X ACIDIFIED	DRAIN	0 K 1 K 2 K 3 K 4 K 0 K 1 K 2 K 3 K 4 K
36.83X BICARB	H ₂ O WATER	36.83X ACIDIFIED	36.83X ACIDIFIED	36.83X ACIDIFIED	DRAIN	0 K 1 K 2 K 3 K 4 K 0 K 1 K 2 K 3 K 4 K
45X BICARB	H ₂ O WATER	45X ACIDIFIED	45X ACIDIFIED	45X ACIDIFIED	DRAIN	2.0 Ca 2.25 Ca 2.5 Ca 2.75 Ca 2.0 Ca 2.25 Ca 2.5 Ca 2.75 Ca
36.1X BICARB	H ₂ O WATER	36.1X ACIDIFIED	36.1X ACIDIFIED	36.1X ACIDIFIED	DRAIN	2.0 Ca 2.25 Ca 2.5 Ca 2.75 Ca


FORM# RPWS-L1

WALL STATION ITEMS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-RPWS-L1	Labels for Recessed Patient Wallstation	1.65
P6FC8	3/8" Tube Compression X 1/2" Male NPT	6.92
P-OR.6	"O"RING for P6FC8 3/8" Tube Compression	0.43
QD-HFCD101235GHT	3/4" MGHT PSF Inline coupling insert w/Shut-off	24.99
QD-HFCD24812	1/2" MNPT PP Inline coupling insert w/Shut-off	22.34
O-RING-WTR	"O" RING for Water Male QD-HFCD24812	0.34
W-ATE2	Wallstation Drain 1" MNPT X 3/4" Barb Elbow	3.70
W-ATA1	Wallstation Drain 1" MNPT X 3/4" Barb Straight	3.10
436-133	Wallstation Drain 1" MNPT x 1 1/2" Socket Schedule 40	4.95
W-M11X1	1 Slot Manifold for Wallstation	10.98
W-M11X2	2 Slot Manifold for Wallstation	12.96
W-M11X3	3 Slot Manifold for Wallstation	15.76
W-M11X4	4 Slot Manifold for Wallstation	26.81
W-M11X6	6 Slot Manifold for Wallstation	28.79

MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com


800.326.5275
 913.438.9700
 FAX 913.438.9701



W-CNC-AGD- AIR GAP DRAIN ASSEMBLY




W-CNC-CAS



AIR GAP DRAIN ASSEMBLY

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-CNC-AGD	Air Gap Drain Assembly	29.95
W-CNC-CAS	Air Gap Body	19.98
W-PL-TAGD	Tee Style Sch 40 Air Gap Drain Assembly	29.95
W-PL-TAGDT	Threaded Tee Style Sch 40 Air Gap Drain Assembly	29.95



W-CAB-282208



WALL STATION CABINET


PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-CAB-282208	Custom Wall Mounting Chase for RPWS (White)	495.25

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com


800.326.5275
913.438.9700
FAX 913.438.9701

CLEAR SCREEN Protection

Protect your new monitor.... Refurbish your old....



Clear Screen for T Machine




CS1002

\$19.95



Fresenius 2008K



B. Braun

"Save your 'T' screen".


"Stops cloudy screens".




medicalsolutionskc.com 800-326-5275

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com


800.326.5275
913.438.9700
FAX 913.438.9701



NOTE COLLAR



SOLENOID
Normally Open



SOLENOID
Normally Closed

BICARB SYSTEM ACCESSORIES

19

SOLENOIDS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
SOL-7524VNC	Solenoid 3/4" 24 VT Normally Closed AC	120.98
SOL-7524VNO	Solenoid 3/4" 24 VT Normally Open AC	193.58
SOL-S201Z1003NC	Stainless Steel 3/4" Normally Closed AC	341.68
SOL-75DC-VNC	Solenoid 3/4" 24 VT Normally Closed DC	120.98
SOL-75DC0VNO	Solenoid 3/4" 24 VT Normally Open DC	193.58
SOL-04-10142	Coil for Zero Differential	26.36
SOL-04-SAN-D	Santilite Diaphragm	32.65
SOL-04-10141	Spring for SS Armature (PLUNGER)	9.85
SOL-04-10296V	SS Amature (PLUNGER)	21.89
SOL-BASE	Black Base under diaphragm 3/4" Thread	89.99
SOL-DENZ	Solenoid Den Connector	6.99
SOL-KIT	Solenoid Santoprene Kit	46.50
SOL-O-RING	"O"RINGS for all Solenoids	5.83

NOTE COLLAR
FOR
NORMALLY
OPEN



*SDS-CB200
2020 or after have
DC solenoids


INJECTORS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
INJ-484	1/2" Injector Venturi-Kynar	121.43



EDUCTORS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ED-4927OPP	Eductor	64.49



SPRAY NOZZLES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
SN-17020	1/4" Spiral Jet Spray Nozzle	27.90
SN-4814T84	1/4" Blue Spray Nozzle 120 degree	27.90



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

BICARB SYSTEM ACCESSORIES

20

HATCH COVER

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
HHC-711	Hinged Hatch Cover-Cam Latch	105.64
HHC-316	Hinged Gasket Cord Kit	13.33
HHC-RETRO KIT	Kit Switching from HHC-912 to HHC-711	46.95



TANK LABELS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
LDT-100	Day Tank Label	11.81
LMT-100	Mix Tank Label	11.81
LDT-50	Day Tank Label	11.81
LMT-50	Mix Tank Label	11.81



KENT METERS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
KM-CA-CONV-KIT	Kent Meter to Carlon Meter	775.00
KM-CA-50	Carlon 1/2" Water Meter	367.33
KM-CA-75	Carlon 3/4" Water Meter	367.33
KM-CA-WL	Wire Lead for Meter	96.25



WAND

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
M-MXR-WD	Wand Assembly for 100 Gallon Tank	137.98
M-MXR-WB	Bicarb Wand Block for 100 Gallon Tank	19.95



HYDROMETERS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
HYDRO-C-08265	Plastic Specific Gravity Hydrometer	165.95
HYDRO-E-6135-6	Hydrometer Jar 500 ml	79.90



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

FUSES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
FU-MXR-20	20 AMP Fuse (120 vt. sys)	39.98
FU-MXR-6.25	6.25 AMP Fuse (24 vt. sys)	39.98
FU-BU-FRN-R-2	Bullet Fuse 20 AMP	28.95
FU-BU-ECNR6.25	Bullet Fuse 6.25 AMP	28.95

SWITCHES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
SWT-24RS	24 Volt Rocker Switch (#1,2,3,4,5)	13.43
SWT-24RS/BR	SMALL 24v Rocker Switch (2020+)	13.43
SWT-120RS	120 Volt Rocker Switch (#6, #7 Switches)	13.43
SWT-120TGL	120 Volt Toggle Switch (not for CB200)	13.43
SWT-1092D5	Green Indicator Light (Chlorine Timer)	13.43
SWT-53413-5	Green Indicator Light (MM1PD Fill Mix)	13.43
SWT-53413-6	Green Indicator Light (MM1PD Recirc, ACID)	13.43
SWT-53415-5	Red Indicator Light (MM1PD Recirc)	13.43
SWT-1092D1	Red Indicator Light (Chlorine Timer)	13.43
SWT-SAFETY	Safety Switch	75.00

SENSORS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
FS-50195k94	Vertical Mount Level Switch	52.46
FS-46515k41	Horizontal Float Switch Side Mount	49.95
FS-135815	Flow Sensor 1"	245.00

PROXIMITY SWITCH

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
SWT-PROX	Proximity Switch-Hartfiel	299.36

CIRCUIT BREAKER

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CB-L881C	Circuit Breaker	38.95



TRANSFORMER

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
TR-50	Transformer 120 VT- 24 AC	65.95
TRC-C-120V	Transformer LVS-150	65.95
SOLA-12549334	SOLA Transformer	195.85



STATIC INLINE MIXER

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
STATIC MIXER	Static Inline Mixer	186.75



ELECTRICAL CONTROL BOXES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
M-RA02	Bicarb Remote Alarm with Reset Mute	641.25
M-MXR-BICFM	Dual Pump Control Box on Single 100 gallon mixers	695.36
PP-24VT	Piezo Buzzer	17.48
M-CT02	Time Delay Chlorine Test Timer	586.95
M-RA03	Remote Alarm for Chlorine Test Timer	129.95
M-RA-W	18 Gauge 2 Conductor Sheathed Wire	.34ft



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

PRESSURE GAUGE

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
PG-B254	Pressure Gauge Brass 0-100 Lower Mount	80.66
PG-AP-0-30	Plastic Pressure Gauge with Subber for APDS-CAB3	261.91
PG-S-5335CBM-100	Pressure Gauge Stainless Steel 0-100 Back Mount	117.49
PG-S-5335CLM-100	Pressure Gauge Stainless Steel 0-100 Lower Mount	117.49
PG-S-254CBM-30	Pressure Gauge Stainless Steel 0-30 Back Mount	117.49
PG-S-254CLM-30	Pressure Gauge Stainless Steel 0-30 Lower Mount	117.49



PRESSURE REGULATOR

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
PR-BHR	Pressure Regulator Bulkhead Nut	9.50
PR-R91S200AA	Pressure Regulator Norgren w/ 1/4" FNPT	63.96
PR-AAP044506	Pressure Regulator Watts w/ 1/4" FNPT	46.95
PR-R362IUY1-S1	Pressure Regulator AOP w/ 1/4" FNPT	55.25



PR-R91S200AA



PR-AAP044506



PR-R362IUY1-S1

CAMLOCKS		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CL-B075	Camlock Female Coupler x 3/4" Male NPT	17.79
CL-B100	Camlock Female Coupler x 1" Male NPT	18.56
CL-C075	Camlock Female Coupler x 3/4" Barb	17.79
CL-C100	Camlock Female Coupler x 1" Barb	18.56
CL-E075	Camlock Male Adapter x 3/4" Barb	10.25
CL-E100	Camlock Male Adapter x 1" Barb	10.58
CL-F075	Camlock Male Adapter x 3/4" Male NPT	10.21
CL-F100	Camlock Male Adapter x 1" Male NPT	10.26
CL-DC075	Camlock Dust Cap 3/4"	17.28
CL-DC100	Camlock Dust Cap 1"	17.77
CL-DP075	Camlock Dust Plug 3/4"	10.33
CL-DP100	Camlock Dust Plug 1"	10.39



Series B



Series C



Series E



Series F



Series DC



Series DP

BULKHEAD FITTINGS		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
BH-50	1/2" TXT Banjo Bulkhead Fitting	27.84
BH-75	3/4" TXT Banjo Bulkhead Fitting	29.62
BH-16439	1/2" White Bulkhead Fitting	39.56



FILTERS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
FIL-CAP 2	.2 Micron Air Filter	116.65
FIL-CAP2/GAC	.2 Micron Air Filter w/Granulated Activated Carbon	127.65
FIL-DCP02000	.2 Micron Air Filter for Bicarb System	149.99
FIL-5-GAC-H	Blue 1/2" Housing for 5" Carbon Filter	29.99
FIL-5-GAC	Granulated Activated Carbon Filter	18.31
FIL-50H-SS	1/2" Clear Bulb, Filter Housing and Filter Screen	38.48
FIL-50-BH	1/2" Clear Bulb	28.50
FIL-50X100	Filter for 1/2" Housing	13.43
FIL-75H-SS	3/4" Clear Bulb, Filter Housing and Filter Screen	65.48
FIL-75-BH	3/4" Clear Bulb	38.50
FIL-75X100	Filter for 3/4" Housing	13.43
FIL-PO.2-10E	.2 Micron 10" RO Storage Tank Air Filter	162.00
FIL-10-H	Housing for 10" 5 micron Air Filters	57.52
FIL-AP-5-10	10" 5 Micron Air Filter	10.55
FIL-AP-5-10/24	10" 5 Micron Air Filter (24 Per Case)	193.13
FIL-AP-5-20	20" 5 Micron Fines Filter (RO Feed)	9.35
FIL-AP-5-20/20	20" 5 Micron Case (20 Per case)	187.06
FIL-PCP-10BB	10" Pentekchlor 20 BB Carbon Block	99.92
FIL-PCP-20BB	20" Pentekchlor 20 BB Carbon Block	219.98
FIL-BB-50233	20" Filter Big Blue Housing	116.51
FIL-BB-AW0296	Wrench SW-3 for BB Housing	9.45
FIL-BB-BB0R	Big Blue filter O Ring	6.35
FIL-GAC20-BB	10" Carbon Block Filter	120.25
TW30-4040	Membranes	495.00



FIL-DCP02000



FIL-5-GAC-H



FIL-5-GAC




FIL-50H-SS



FIL-75H-SS

IWAKI PUMP		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
PMP-IW-55RLT	Iwaki Pump 55 RLT (complete motor and wet end kit)	795.75
PMP-IW-55RLT-WEK	Iwaki Pump Wet End Kit 55 RLT	517.20
PMP-IW-70RLZT	Iwaki Pump 70 RLZT	981.76
PMP-IW-70RLZT-WEK	Iwaki Pump Wet End Kit-70 RLZT	586.58
PMP-IW-70RLT	Iwaki Pump 70 RLT	981.76
PMP-IW-70RLT-WEK	Iwaki Pump Wet End Kit 70 RLT	586.58
PMP-IW-30RLT	Iwaki Pump 30 RLT	524.43
PMP-IW-100RLT	Iwaki Pump 100 RLT	798.75




MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

MARCH PUMP TE-7R

PART NUMBER	DIAGRAM #	PRODUCT DESCRIPTION	PRICE \$
PMP-M-TE7R		TE7R Large March Pump (Motor and Wet End Kit)	1199.00
PMP-M-TE7R-WEK	#1- #8	TE7R Wet End Kit	576.30
PMP-M-0011-10	#3	Front Housing TE7R	228.25
PMP-M-0010-10	#4	Viton O-RING TE7R	33.75
PMP-M-0009-10	#5	Ceramic Washer TE7R	12.19
PMP-M-0039-10	#6	Ceramic Spindle TE7R	49.98
PMP-M-0159-05	#7	Impellar Assembly TE7R	289.58
PMP-M-0067-01	#8	Rear Magnet Housing TE7R	193.55
PMP-M-0124-01	ALT #8	Rear Magnet Housing Kynar	215.55
PMP-M-0130-02	#9	Drive Magnet Housing Assembly TE7R	429.17
PMP-M-0092-01	#10	Motor to Pump Bracket TE7R	141.75
PMP-M-0016-14	#11	Motor TE7R-DOES NOT INCLUDE WET END KIT	883.44



PUMPS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
PMP-M-AC4C	March Pump AC4C-MD	467.50
PMP-SF-2088	Shur Flo Pump 3.3 GPM 115VT	161.84
PMP-SF-2088-PA	Shur Flo Pump including purge assembly	391.15



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

HOSE AND HOSE CLAMPS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
HC-100	Hose Clamp for 1" Hose OD	2.74
HC-3/4	Hose Clamp for 3/4" Hose OD	2.74
HC-1/2	Hose Clamp for 1/2" Hose OD	2.74
HC-3/8	Hose Clamp for 3/8" Hose OD	2.74
HC-BLACK	Black Hose Clamp for 3/8" OD	2.38
HC-BLUE	Blue Hose Clamp for 3/8" OD	2.38
H-PB-1-G	1" Poly Braided Hose Gray	2.74
H-PB-3/4-G	3/4" Poly Braided Hose Gray	4.34
H-PB-1/2-G	1/2" Poly Braided Hose Gray	2.74
H-PB-1/2-CL	1/2" Poly Braided Hose Clear	2.98
H-PB-3/8-G	3/8" Poly Braided Hose Gray	2.74
H-PB-3/8-CL	3/8" Poly Braided Hose Clear	2.98
H-PV-3/4-CL	3/4" Clear Poly Vinyl Hose	2.90
H-PV-3/8	3/8" Clear Poly Vinyl Hose	2.74



CLICKS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CLIC-12	#12 Click for Drain Hose	1.85
CLIC-15	#15 Click for Drain Hose	1.85
CLIC-17	#17 Click for 3/8" PVC	1.95
CLIC-20	#20 Click for 1/2" PVC	2.78
CLIC-25	#25 Click for 3/4" PVC	3.44
CLIC-32	#32 Click for 1" PVC	4.86
CLIC-47	#47 Click for 1-1/2" PVC	7.76



TUBING		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
1/2" BLUE	1/2" Poly Tubing Blue	0.66
1/2" ORANGE	1/2" Poly Tubing Orange (not stock)	0.66
1/2" RED	1/2" Poly Tubing Red	0.66
1/2" YELLOW	1/2" Poly Tubing Yellow	0.66
1/4" BLACK	1/4" Poly Tubing Black	0.66
1/4" BLUE	1/4" Poly Tubing Blue	0.66
1/4" RED	1/4" Poly Tubing Red	0.66
1/4" YELLOW	1/4" Poly Tubing Yellow	0.66
1/4" ORANGE	1/4" Poly Tubing Orange	0.66
1/4" WHITE	1/4" Poly Tubing White	0.66
3/8" BLUE	3/8" Poly Tubing Blue	0.65
3/8" GREEN	3/8" Poly Tubing Green	0.65
3/8" ORANGE	3/8" Poly Tubing Orange	0.65
3/8" RED	3/8" Poly Tubing Red	0.65
3/8" WHITE	3/8" Poly Tubing White	0.65
3/8" YELLOW	3/8" Poly Tubing Yellow	0.65
5/8" YELLOW	5/8" Poly Tubing Yellow	0.85
5/8" BLUE	5/8" Poly Tubing Blue	0.85
5/8" RED	5/8" Poly Tubing Red	0.85



COILED CONCENTRATED TUBING		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
C80550	Red Coiled Concentrated Tubing	16.95
C80553	Blue Coiled Concentrated Tubing	16.95
C80558	Clear Coiled Concentrated Tubing	16.95



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

COUPLING SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
829-002	1/4" SXS Coupling	8.96
829-003	3/8" SXS Coupling	8.94
829-005	1/2" SXS Coupling	5.87
829-007	3/4" SXS Coupling	7.68
829-010	1" SXS Coupling	7.95
829-012	1-1/4" SXS Coupling	12.45
829-015	1-1/2" SXS Coupling	13.92
829-020	2" SXS Coupling	14.73

**FEMALE ADAPTER SOC X FIPT**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
835-002	1/4" SXT Female Adapter	9.95
835-003	3/8" SXT Female Adapter	10.26
835-005	1/2" SXT Female Adapter	6.45
835-007	3/4" SXT Female Adapter	9.98
835-010	1" SXT Female Adapter	13.79
835-012	1-1/4" SXT Female Adapter	17.74
835-015	1-1/2" SXT Female Adapter	27.56
835-020	2" SXT Female Adapter	45.26

**MALE ADAPTER MIPT X SOC**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
836-005	1/2" SXT Male Adapter	7.96
836-007	3/4" SXT Male Adapter	7.57
836-010	1" SXT Male Adapter	13.98
836-012	1-1/4" SXT Male Adapter	17.74
836-015	1-1/2" SXT Male Adapter	19.96
836-020	2" SXT Male Adapter	28.94



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203

www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

COUPLING FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
830-002	1/4" TXT Coupling	6.29
830-003	3/8" TXT Coupling	6.56
830-005	1/2" TXT Coupling	7.58
830-007	3/4" TXT Coupling	11.13
830-010	1" TXT Coupling	11.86
830-012	1-1/4" TXT Coupling	14.96
830-015	1-1/2" TXT Coupling	35.68
830-020	2" TXT Coupling	32.89




90° ELBOW SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
806-002	1/4" SXS 90° Elbow	4.98
806-003	3/8" SXS 90° Elbow	4.98
806-005	1/2" SXS 90° Elbow	3.94
806-007	3/4" SXS 90° Elbow	3.97
806-010	1" SXS 90° Elbow	6.18
806-012	1-1/4" SXS 90° Elbow	8.71
806-015	1-1/2" SXS 90° Elbow	9.58
806-020	2" SXS 90° Elbow	15.53



90° ELBOW LONG SWEEP SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
806-005LSF	1/2" SXS 90° Long Sweep	25.51
806-007LSF	3/4" SXS 90° Long Sweep	28.63
806-010LSF	1" SXS 90° Long Sweep	32.92
806-012LSF	1-1/4" SXS 90° Long Sweep	54.13
806-015LSF	1-3/4" SXS 90° Long Sweep	54.82
806-020LSF	2" SXS 90° Long Sweep	59.27



90° ELBOW SWEEP SOC X SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
806-005S	1/2" SXS 90° Sweep	7.36
806-007S	3/4" SXS 90° Sweep	14.39
806-010S	1" SXS 90° Sweep	15.35
806-012S	1-1/4" SXS 90° Sweep	19.93
806-015S	1-3/4" SXS 90° Sweep	21.82
806-020S	2" SXS 90° Sweep	26.43



90° SOC X FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
807-002	1/4" SXT 90° Elbow	12.93
807-003	3/8" SXT 90° Elbow	11.61
807-005	1/2" SXT 90° Elbow	12.77
807-007	3/4" SXT 90° Elbow	14.99
807-010	1" SXT 90° Elbow	16.38
807-012	1-1/4" SXT 90° Elbow	19.40
807-015	1-1/2" SXT 90° Elbow	23.26
807-020	2" SXT 90° Elbow	34.99



90° ELBOW FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
808-002	1/4" TXT 90° Elbow	6.98
808-003	3/8" TXT 90° Elbow	7.32
808-005	1/2" TXT 90° Elbow	7.51
808-007	3/4" TXT 90° Elbow	10.89
808-010	1" TXT 90° Elbow	11.75
808-012	1-1/4" TXT 90° Elbow	13.56
808-015	1-1/2" TXT 90° Elbow	18.75
808-020	2" TXT 90° Elbow	22.56



90° STREET ELBOW SPRIG X SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
809-005	1/2" 90° Street Elbow	4.96
809-007	3/4" 90° Street Elbow	5.87
809-010	1" 90° Street Elbow	9.25
809-012	1-1/4" 90° Street Elbow	12.91
809-015	1-1/2" 90° Street Elbow	11.73
809-020	2" 90° Street Elbow	15.88

**45° ELBOW SOC**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
817-002	1/4" SXS 45° Elbow	11.43
817-003	3/8" SXS 45° Elbow	11.17
817-005	1/2" SXS 45° Elbow	5.56
817-007	3/4" SXS 45° Elbow	8.94
817-010	1" SXS 45° Elbow	14.98
817-012	1-1/4" SXS 45° Elbow	19.34
817-015	1-1/2" SXS 45° Elbow	23.49
817-020	2" SXS 45° Elbow	29.94

**45° ELBOW FIPT**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
819-002	1/4" TXT 45° Elbow	15.57
819-003	3/8" TXT 45° Elbow	16.92
819-005	1/2" TXT 45° Elbow	15.96
819-007	3/4" TXT 45° Elbow	16.85
819-010	1" TXT 45° Elbow	17.57
819-012	1-1/4" TXT 45° Elbow	17.92
819-015	1-1/2" TXT 45° Elbow	27.46
819-020	2" TXT 45° Elbow	42.20



45° STREET ELBOW SPIG X SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
827-005	1/2" 45° Street Elbow	7.56
827-007	3/4" 45° Street Elbow	11.49
827-010	1" 45° Street Elbow	14.73
827-012	1-1/4" 45° Street Elbow	17.77
827-015	1-1/2" 45° Street Elbow	20.96
827-020	2" 45° Street Elbow	25.85

**TEE SOC**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
801-002	1/4" SXSXS Tee	9.96
801-003	3/8" SXSXS Tee	9.96
801-005	1/2" SXSXS Tee	9.35
801-007	3/4" SXSXS Tee	9.56
801-010	1" SXSXS Tee	11.98
801-012	1-1/4" SXSXS Tee	31.96
801-015	1-1/2" SXSXS Tee	31.96
801-020	2" SXSXS Tee	37.78

**TEE SOC X FIPT**

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
802-002	1/4" SXSXT Tee	18.24
802-005	1/2" SXSXT Tee	14.26
802-007	3/4" SXSXT Tee	18.63
802-010	1" SXSXT Tee	21.93
802-012	1-1/4" SXSXT Tee	21.31
802-015	1-1/2" SXSXT Tee	25.28
802-020	2" SXSXT Tee	30.87



TEE FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
805-002	1/4" TXTXT Tee	10.58
805-003	3/8" TXTXT Tee	13.45
805-005	1/2" TXTXT Tee	10.58
805-007	3/4" TXTXT Tee	11.78
805-010	1" TXTXT Tee	24.52
805-012	1-1/4" TXTXT Tee	29.38
805-015	1-1/2" TXTXT Tee	29.84
805-020	2" TXTXT Tee	34.95



CROSS SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
820-002	1/4" SXS Cross	44.92
820-005	1/2" SXS Cross	44.98
820-007	3/4" SXS Cross	43.17
820-010	1" SXS Cross	57.56
820-012	1-1/4" SXS Cross	58.56
820-015	1-1/2" SXS Cross	65.59
820-020	2" SXS Cross	79.97



CAP SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
847-002	1/4" Cap Socket	8.79
847-003	3/8" Cap Socket	8.79
847-005	1/2" Cap Socket	7.23
847-007	3/4" Cap Socket	7.51
847-010	1" Cap Socket	11.79
847-012	1-1/4" Cap Socket	13.81
847-015	1-1/2" Cap Socket	13.81
847-020	2" Cap Socket	23.31



PVC SCHEDULE 80 FITTINGS

36

CAP FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
848-002	1/4" Threaded Cap	7.53
848-003	3/8" Threaded Cap	7.53
848-005	1/2" Threaded Cap	9.38
848-007	3/4" Threaded Cap	10.57
848-010	1" Threaded Cap	12.56
848-012	1-1/4" Threaded Cap	12.98
848-015	1-1/2" Threaded Cap	14.54
848-020	2" Threaded Cap	19.64



PLUG SPIG

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
849-005	1/2" SPIG Plug	6.89
849-007	3/4" SPIG Plug	7.44
849-010	1" SPIG Plug	8.14
849-012	1-1/4" SPIG Plug	9.83
849-015	1-1/2" SPIG Plug	10.59
849-020	2" SPIG Plug	13.36



PLUG MIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
850-002	1/4" Threaded Plug	8.46
850-003	3/8" Threaded Plug	7.84
850-005	1/2" Threaded Plug	7.46
850-007	3/4" Threaded Plug	7.75
850-010	1" Threaded Plug	8.42
850-012	1-1/4" Threaded Plug	13.81
850-015	1-1/2" Threaded Plug	15.53
850-020	2" Threaded Plug	14.76



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

REDUCER BUSHING FLUSH SPIG X SOC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
RB-837-052	3/8" X 1/4" SXS Reducer Bushing	5.72
RB-837-072	1/2" X 1/4" SXS Reducer Bushing	5.72
RB-837-073	1/2" X 3/8" SXS Reducer Bushing	5.72
RB-837-101	3/4" X 1/2" SXS Reducer Bushing	4.89
RB-837-129	1" X 3/8" SXS Reducer Bushing	9.83
RB-837-130	1" X 1/2" SXS Reducer Bushing	5.34
RB-837-131	1" X 3/4" SXS Reducer Bushing	5.42
RB-837-166	1-1/4" X 1/2" SXS Reducer Bushing	7.73
RB-837-167	1-1/4" X 3/4" SXS Reducer Bushing	7.73
RB-837-168	1-1/4" X 1" SXS Reducer Bushing	8.21
RB-837-209	1-1/2" X 1/2" SXS Reducer Bushing	10.53
RB-837-210	1-1/2" X 3/4" SXS Reducer Bushing	10.53
RB-837-211	1-1/2" X 1" SXS Reducer Bushing	11.58
RB-837-212	1-1/2" X 1-1/4" SXS Reducer Bushing	10.53
RB-837-247	2" X 1/2" SXS Reducer Bushing	14.93
RB-837-248	2" X 3/4" SXS Reducer Bushing	14.93
RB-837-249	2" X 1" SXS Reducer Bushing	14.93
RB-837-250	2" X 1-1/4" SXS Reducer Bushing	14.93
RB-837-251	2" X 1/2" SXS Reducer Bushing	14.93
RB-837-290	2-1/2" X 1-1/4" SXS Reducer Bushing	25.95
RB-837-291	2-1/2" X 1-1/2" SXS Reducer Bushing	25.95
RB-837-292	2-1/2" X 2" SXS Reducer Bushing	25.95



REDUCER BUSHING FLUSH SPIG X FIPT		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
RB-838-052	3/8"X 1/4" SXT Reducer Bushing	4.72
RB-838-071	1/2"X 1/8" SXT Reducer Bushing	4.72
RB-838-072	1/2"X 1/4" SXT Reducer Bushing	5.79
RB-838-073	1/2"X 3/8" SXT Reducer Bushing	5.72
RB-838-098	3/4"X 1/4" SXT Reducer Bushing	4.56
RB-838-099	3/4"X 3/8" SXT Reducer Bushing	7.93
RB-838-101	3/4"X 1/2" SXT Reducer Bushing	4.56
RB-838-127	1"X 1/8" SXT Reducer Bushing	4.93
RB-838-129	1"X 3/8" SXT Reducer Bushing	9.85
RB-838-130	1"X 1/2" SXT Reducer Bushing	5.32
RB-838-131	1"X 3/4" SXT Reducer Bushing	5.42
RB-838-166	1-1/4"X 1/2" SXT Reducer Bushing	9.92
RB-838-167	1-1/4"X 3/4" SXT Reducer Bushing	9.92
RB-838-168	1-1/4"X 1" SXT Reducer Bushing	9.92
RB-838-206	1-1/2"X 1/8" SXT Reducer Bushing	12.54
RB-838-207	1-1/2"X 1/4" SXT Reducer Bushing	12.91
RB-838-209	1-1/2"X 1/2" SXT Reducer Bushing	15.14
RB-838-210	1-1/2"X 3/4" SXT Reducer Bushing	12.54
RB-838-211	1-1/2"X 1" SXT Reducer Bushing	13.79
RB-838-212	1-1/2"X 1-1/4" SXT Reducer Bushing	12.54
RB-838-247	2"X 1/2" SXT Reducer Bushing	15.34
RB-838-248	2"X 3/4" SXT Reducer Bushing	16.18
RB-838-249	2"X 1" SXT Reducer Bushing	18.18
RB-838-250	2"X 1-1/4" SXT Reducer Bushing	18.34
RB-838-251	2"X 1-1/2" SXT Reducer Bushing	18.34





REDUCER BUSHING FLUSH MIPT X FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
TRB-839-041	1/4" X 1/8" TXT Reducer Bushing	11.84
TRB-839-052	3/8" X 1/4" TXT Reducer Bushing	12.10
TRB-839-072	1/2" X 1/4" TXT Reducer Bushing	11.84
TRB-839-073	1/2" X 3/8" TXT Reducer Bushing	12.18
TRB-839-098	3/4" X 1/4" TXT Reducer Bushing	8.77
TRB-839-099	3/4" X 3/8" TXT Reducer Bushing	9.30
TRB-839-101	3/4" X 1/2" TXT Reducer Bushing	6.61
TRB-839-128	1" X 1/4" TXT Reducer Bushing	14.23
TRB-839-129	1" X 3/8" TXT Reducer Bushing	15.97
TRB-839-130	1" X 1/2" TXT Reducer Bushing	11.39
TRB-839-131	1" X 3/4" TXT Reducer Bushing	11.39
TRB-839-166	1-1/4" X 1/2" TXT Reducer Bushing	12.32
TRB-839-167	1-1/4" X 3/4" TXT Reducer Bushing	11.98
TRB-839-168	1-1/4" X 1" TXT Reducer Bushing	12.94
TRB-839-209	1-1/2" X 1/2" TXT Reducer Bushing	15.96
TRB-839-210	1-1/2" X 3/4" TXT Reducer Bushing	14.26
TRB-839-211	1-1/2" X 1" TXT Reducer Bushing	14.26
TRB-839-212	1-1/2" X 1-1/4" TXT Reducer Bushing	14.26
TRB-839-247	2" X 1/2" TXT Reducer Bushing	18.98
TRB-839-248	2" X 3/4" TXT Reducer Bushing	19.98
TRB-839-249	2" X 1" TXT Reducer Bushing	18.98
TRB-839-250	2" X 1-1/4" TXT Reducer Bushing	18.98
TRB-839-251	2" X 1-1/2" TXT Reducer Bushing	18.98
TRB-839-292	2-1/2" X 2" TXT Reducer Bushing	25.81



UNION (OLD STYLE) FIPT W/EDPM O-RING SEAL

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
898-005	1/2" Threaded Union EPDM	16.68
898-007	3/4" Threaded Union EPDM	21.88
898-010	1" Threaded Union EPDM	28.70
898-012	1-1/4" Socket Union EDPM	55.99
898-015	1-1/2" Socket Union EDPM	58.69
898-020	2" Socket Union EDPM	93.89




UNION (OLD STYLE) SOC W/EDPM O-RING SEAL

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
897-005	1/2" Socket Union EPDM	13.57
897-007	3/4" Socket Union EPDM	17.23
897-010	1" Socket Union EPDM	19.99
897-012	1-1/4" Socket Union EDPM	39.14
897-015	1-1/2" Socket Union EDPM	44.40
897-020	2" Socket Union EDPM	47.84
897-025	2-1/2" Socket Union EDPM	76.05
897-030	3" Socket Union EDPM	89.02
897-040	4" Socket Union EDPM	123.13




TANK ADAPTER WITH EDPM GASKET SOC X F

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
871-E-005	1/2" SXT Bulkhead	54.79
871-E-007	3/4" SXT Bulkhead	55.65
871-E-010	1" SXT Bulkhead	52.99
871-E-012	1-1/4" SXT Bulkhead	70.88
871-E-015	1-1/2" SXT Bulkhead	75.66
871-E-020	2" SXT Bulkhead	91.91



TANK ADAPTER WITH EDPM GASKET FIPT

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
872-E-005	1/2" TXT Bulkhead	62.42
872-E-007	3/4" TXT Bulkhead	66.05
872-E-010	1" TXT Bulkhead	74.71
872-E-012	1-1/4" TXT Bulkhead	71.91
872-E-015	1-1/2" TXT Bulkhead	77.78
872-E-020	2" TXT Bulkhead	101.74



BULKHEAD FITTINGS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
BH-50	1/2" TXT Banjo Bulkhead Fitting	27.84
BH-75	3/4" TXT Banjo Bulkhead Fitting	29.62
BH-16439	1/2" White Bulkhead	39.56



NIPPLES 1/4"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
880-005	Close 1/4" NIPPLE	3.49
880-015	1-1/2" Long 1/4" Nipple	4.33
880-020	2" Long 1/4" Nipple	4.33
880-025	2-1/2" Long 1/4" Nipple	5.75
880-030	3" Long 1/4" Nipple	5.75



NIPPLES 3/8"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
881-005	Close 3/8" Nipple	5.88
881-015	1-1/2" Long 3/8" Nipple	5.92
881-020	2" Long 3/8" Nipple	5.98
881-025	2-1/2" Long 3/8" Nipple	6.47
881-030	3" Long 3/8" Nipple	6.47



NIPPLES 1/2"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
882-005	Close 1/2" Nipple	3.14
882-015	1-1/2" Long 1/2" Nipple	3.15
882-020	2" Long 1/2" Nipple	3.41
882-025	2-1/2" Long 1/2" Nipple	3.98
882-030	3" Long 1/2" Nipple	3.76
882-035	3-1/2" Long 1/2" Nipple	3.76
882-040	4" Long 1/2" Nipple	4.47
882-045	4-1/2" Long 1/2" Nipple	4.47
882-050	5" Long 1/2" Nipple	4.55
882-055	5-1/2" Long 1/2" Nipple	4.47
882-060	6" Long 1/2" Nipple	5.50
882-060M	6" Long 1/2" Nipple Machine Cut	6.35



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

NIPPLES 3/4"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
883-005	Close 3/4" Nipple	2.99
883-015	1-1/2" Long 3/4" Nipple	3.48
883-020	2" Long 3/4" Nipple	3.40
883-025	2-1/2" Long 3/4" Nipple	3.75
883-030	3" Long 3/4" Nipple	3.75
883-035	3-1/2" Long 3/4" Nipple	5.15
883-040	4" Long 3/4" Nipple	5.18
883-045	4-1/2" Long 3/4" Nipple	5.63
883-050	5" Long 3/4" Nipple	6.23
883-055	5-1/2" Long 3/4" Nipple	6.45
883-060	6" Long 3/4" Nipple	6.53
883-140M	14" Long x 3/4" Nipple Machine C	9.72



NIPPLES 1"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
884-005	Close 1" Nipple	4.49
884-015	1-1/2" Long 1" Nipple	6.33
884-020	2" Long 1" Nipple	6.33
884-025	2-1/2" Long 1" Nipple	6.40
884-030	3" Long 1" Nipple	6.40
884-035	3-1/2" Long 1" Nipple	6.54
884-040	4" Long 1" Nipple	6.54
884-045	4-1/2" Long 1" Nipple	8.89
884-050	5" Long 1" Nipple	8.89
884-055	5-1/2" Long 1" Nipple	9.72
884-060	6" Long 1" Nipple	9.99



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

NIPPLES 1-1/4"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
885-005	Close 1-1/4" Nipple	6.52
885-015	1-1/2" Long 1-1/4" Nipple	5.34
885-020	2" Long 1-1/4" Nipple	6.75
885-025	2-1/2" Long 1-1/4" Nipple	7.39
885-030	3" Long 1-1/4" Nipple	7.53
885-035	3-1/2" Long 1-1/4" Nipple	7.72
885-040	4" Long 1-1/4" Nipple	7.72
885-045	4-1/2" Long 1-1/4" Nipple	7.97
885-050	5" Long 1-1/4" Nipple	7.97
885-055	5-1/2" Long 1-1/4" Nipple	8.48
885-060	6" Long 1-1/4" Nipple	8.48



NIPPLES 1-1/2"

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
886-005	Close 1-1/2" Nipple	6.99
886-015	1-1/2" Long 1-1/2" Nipple	6.99
886-020	2" Long 1-1/2" Nipple	7.80
886-025	2-1/2" Long 1-1/2" Nipple	8.98
886-030	3" Long 1-1/2" Nipple	8.98
886-035	3-1/2" Long 1-1/2" Nipple	10.42
886-040	4" Long 1-1/2" Nipple	11.93
886-045	4-1/2" Long 1-1/2" Nipple	13.55
886-050	5" Long 1-1/2" Nipple	13.91
886-055	5-1/2" Long 1-1/2" Nipple	12.43
886-060	6" Long 1-1/2" Nipple	15.34



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

BALL VALVES		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
BV-A1610-015	1-1/2" ASHAI True Union Ball Valve	140.20
BV-6621-007	3/4" Compact Ball Valve TXT	69.22
BV-2121-010	1" Compact Ball Valve TXT	75.72
BV-6622-007	3/4" Compact Ball Valve SXS	69.22
BV-2122-010	1" Compact Ball Valve SXS	75.72
BV-2122-015	1-1/2" Compact Ball Valve SXS	82.35
BV-2264F4MB	1/4" Female Thread x Male Barb Labcock Ball Valve	51.23
BV-2621-005G	1/2" Utility Ball Valve TXT (GRAY)	28.56
BV-2621-005W	1/2" Sample Port Ball Valve TXT (WHITE)	28.56
BV-2621-SAMPLE	1/2" Sample Port Ball Valve Set up	35.26
BV-ASHAI-3WAY	ASHAI 3 WAY 1/2" Ball Valve Tee TYPE 23	188.12
BV-PVC-350E	Ball Valve 1/2" 3 WAY (w/gray handle)	57.64
BV-SMC-25FF	Ball Valve 1/4" Female NPT x 1/4" Female NPT	22.99
BV-SMC-25MF	Ball Valve 1/4" Male NPT x 1/4" Female NPT	22.99
BV-SS-4F	Ball Valve Stainless Steel 1/4" Male x Female NPT	23.50



CHECK VALVES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
CV-4529-005	1/2" Check Valve True Union Ball Check	111.58
CV-4529-007	3/4" Check Valve True Union Ball Check	126.40
CV-4529-010	1" Check Valve True Union Ball Check	117.14
CV-4529-015	1-1/2" Check Valve True Union Ball Check	174.99
CV-6858F8F-E	1/2" Check Valve Black	22.85
CV-68712F12F	3/4" Check Valve Black	25.95



GATE VALVES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
GV-2021-005	1/2" Threaded Gate Valve	59.99
GV-2021-010	1" Threaded Gate Valve	98.53
GV-2022-005	1/2" Socket Gate Valve	58.26
GV-2022-010	1" Socket Gate Valve	98.53



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

BARBED FITTINGS

47

PVC

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
HB-P4MCB4	1/4" Male NPT x 1/4" Barb (black)	1.48
HB-P4MEB4	1/4" Male NPT x 1/4" Barb Elbow (black)	1.94
1407-010	1" Female NPT X 1" Barb Elbow	4.89
1407-131	3/4" Male NPT X 1" Barb Elbow	5.39
1413-007	3/4" Male NPT X 3/4" Barb Elbow	5.56
1413-010	1" Male NPT X 1" Barb Elbow	5.71
1413-101	1/2" Male NPT X 3/4" Barb Elbow	6.72
1435-007	3/4" Female NPT X 3/4" Barb Elbow	5.53
1436-007	3/4" Male NPT X 3/4" Barb Straight	7.82
1436-010	1" Male NPT X 1" Barb Straight	8.73
1436-074	1/2 Male NPT X 3/4" Barb Straight	5.78
1436-101	3/4" Male NPT x 1/2" Barb Straight	5.35
1436-102	3/4" Male NPT X 1" Barb Straight	4.96
1436-131	1" Male NPT X 3/4" Barb Straight	4.89



POLYPROPYLENE FITTINGS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
POLY-64822	1/2" Male NPT X 3/8" Barb Straight	1.82
POLY-64489	1/2" Male NPT X 3/8" Barb Elbow	1.82

CUSTOM FAUCET

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
PVC-80-005	1/2" Formed Faucet	14.77



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701



GARDEN HOSE THREAD BODY POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD101235GHT	3/8"	3/4" MGHT PSF Inline coupling insert w/Shut-off	24.99
QD-HFCD191235GHT	3/8"	3/4" FGHT PSF Inline coupling insert w/Shut-off	28.68
QD-HFC101235GHT	3/8"	3/8" MGHT PSF Coupler body Straight-thru	25.74
QD-HFC191235GHT	3/8"	3/4" FGHT PSF Coupler body Straight-thru	24.94



HIGH FLOW COUPLINGS BODIES POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD10612	3/8"	3/8" MNPT PP Coupler body w/Shut-off	29.99
QD-HFCD10812	3/8"	1/2" MNPT PP Coupler body w/Shut-off	34.95
QD-HFC10612	3/8"	3/8" MNPT PSF Coupler body Straight-thru	25.96
QD-HFC10812	3/8"	1/2" MNPT PSF Coupler body Straight-thru	24.50



HIGH FLOW COUPLINGS BODIES POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD10635	3/8"	3/8" MNPT PSF Coupler body w/Shut-off	26.93
QD-HFCD10835	3/8"	1/2" MNPT PSF Coupler body w/Shut-off	34.95
QD-HFC10635	3/8"	3/8" MNPT PSF Coupler body Straight-thru	35.95
QD-HFC10835	3/8"	1/2" MNPT PSF Coupler body Straight-thru	35.22



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

QUICK DISCONNECT FITTINGS

49

IN-LINE HOSE BARB BODY POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD17612	3/8"	3/8" ID PP Barbed Inline body w/Shut-off	32.36
QD-HFCD17812	3/8"	1/2" ID PP Barbed Inline body w/Shut-off	32.36
QD-HFCD171212	3/8"	3/4" ID PP Barbed Inline body w/Shut-off	32.85
QD-HFC17612	3/8"	3/8" ID PP Barbed inline body Straight-thru	26.36
QD-HFC17812	3/8"	1/2" ID PP Barbed inline body Straight-thru	26.95
QD-HFC171212	3/8"	3/4" ID PP Barbed inline body Straight-thru	28.12



IN-LINE HOSE BARB BODY POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD17635	3/8"	3/8" ID PSF Barbed Inline body w/Shut-off	38.42
QD-HFCD17835	3/8"	1/2" ID PSF Barbed Inline body w/Shut-off	38.42
QD-HFCD171235	3/8"	3/4" ID PSF Barbed Inline body w/Shut-off	38.42



PANEL MOUNT HOSE BARB BODY POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD16612NA	3/8"	3/8" ID PP Panel mount body w/Shut-off	34.95
QD-HFCD16812NA	3/8"	1/2" ID PP Panel mount body w/Shut-off	38.42
QD-HFCD161212NA	3/8"	3/4" ID PP Panel mount body w/Shut-off	34.95
QD-HFC16612NA	3/8"	3/8" ID PP Panel mount body Straight-thru	29.95
QD-HFC16812NA	3/8"	1/2" ID PP Panel mount body Straight-thru	29.95
QD-HFC161212NA	3/8"	3/4" ID PP Panel mount body Straight-thru	29.95



PANEL MOUNT HOSE BARB BODY POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD16635NA	3/8"	3/8" ID PSF Panel mount body w/Shut-off	38.96
QD-HFCD16835NA	3/8"	1/2" ID PSF Panel mount body w/Shut-off	38.42
QD-HFCD161235NA	3/8"	3/4" ID PSF Panel mount body w/Shut-off	38.42
QD-HFC16635NA	3/8"	3/8" ID PSF Panel mount body Straight-thru	39.76
QD-HFC16835NA	3/8"	1/2" ID PSF Panel mount body Straight-thru	39.88
QD-HFC161235NA	3/8"	3/4" ID PSF Panel mount body Straight-thru	38.65



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

QUICK DISCONNECT FITTINGS

50

HIGH FLOW COUPLING INSERTS

GARDEN HOSE THREAD INSERT POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD261235GHT	3/8"	3/4" FGHT PSF Inline coupling insert w/Shut-off	32.26
QD-HFC241235GHT	3/8"	3/4" MGHT PSF Inline coupling insert Straight-thru	17.26
QD-HFC261235GHT	3/8"	3/4" FGHT PSF Inline coupling insert Straight-thru	24.92
QD-HFCD241235GHT	3/8"	3/4" MGHT PSF Inline coupling insert w/Shut-off	32.26



PIPE THREAD INSERT POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD24612	3/8"	3/8" MNPT PP Inline coupling insert w/Shut-off	25.18
QD-HFCD24812	3/8"	1/2" MNPT PP Inline coupling insert w/Shut-off	22.24
QD-HFC24612	3/8"	3/8" MNPT PP Inline coupling insert Straight-thru	10.92
QD-HFC24812	3/8"	1/2" MNPT PP Inline coupling insert Straight-thru	10.92



PIPE THREAD INSERT POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD24635	3/8"	3/8" MNPT PSF Inline coupling insert w/Shut-off	29.95
QD-HFCD24835	3/8"	1/2" MNPT PSF Inline coupling insert w/Shut-off	29.95
QD-HFC24635	3/8"	3/8" MNPT PSF Inline coupling insert Straight-thru	19.68
QD-HFC24835	3/8"	1/2" MNPT PSF Inline coupling insert Straight-thru	19.68



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

QUICK DISCONNECT FITTINGS

51

HIGH FLOW COUPLING INSERTS

HOSE BARB INSERT POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD22612	3/8"	3/8" ID PP Barbed coupling insert w/Shut-off	24.72
QD-HFCD22812	3/8"	1/2" ID PP Barbed coupling insert w/Shut-off	24.95
QD-HFCD221212	3/8"	3/4" ID PP Barbed coupling insert w/Shut-off	25.50
QD-HFC22612	3/8"	3/8" ID PP Barbed coupling insert Straight-thru	11.52
QD-HFC22812	3/8"	1/2" ID PP Barbed coupling insert Straight-thru	16.56
QD-HFC221212	3/8"	3/4" ID PP Barbed coupling insert Straight-thru	12.80



HOSE BARB INSERT POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD22635	3/8"	3/8" ID PSF Barbed coupling insert w/Shut-off	28.22
QD-HFCD22835	3/8"	1/2" ID PSF Barbed coupling insert w/Shut-off	28.85
QD-HFCD221235	3/8"	3/4" ID PSF Barbed coupling insert w/Shut-off	29.46
QD-HFC22635	3/8"	3/8" ID PSF Barbed coupling insert Straight-thru	10.52
QD-HFC22835	3/8"	1/2" ID PSF Barbed coupling insert Straight-thru	16.62
QD-HFC221235	3/8"	3/4" ID PSF Barbed coupling insert Straight-thru	18.76



ELBOW HOSE BARB INSERT POLYPROPYLENE (gray)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD23612	3/8"	3/8" ID PP Barbed coupling insert w/Shut-off	29.23
QD-HFCD23812	3/8"	1/2" ID PP Barbed coupling insert w/Shut-off	28.78
QD-HFC23612	3/8"	3/8" ID PP Barbed elbow insert Straight-thru	18.47
QD-HFC23812	3/8"	1/2" ID PP Barbed elbow insert Straight-thru	18.78



ELBOW HOSE BARB INSERT POLYSULFONE (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-HFCD23635	3/8"	3/8" ID PSF Barbed coupling insert w/Shut-off	35.22
QD-HFCD23835	3/8"	1/2" ID PSF Barbed coupling insert w/Shut-off	31.56
QD-HFC23635	3/8"	3/8" ID PSF Barbed elbow insert Straight-thru	19.44
QD-HFC23835	3/8"	1/2" ID PSF Barbed elbow insert Straight-thru	19.35



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

QUICK DISCONNECT FITTINGS

52

LOW FLOW COUPLING INSERTS

MALE PIPE THREAD ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC1002	1/8"	1/8" MT Pipe thread body Straight-thru	17.42
QD-PMC1004	1/8"	1/4" MT Pipe thread body Straight-thru	13.55
QD-PLC10004	1/4"	1/4" MT Pipe thread body Straight-thru	25.50
QD-PLC10006	1/4"	3/8" MT Pipe thread body Straight-thru	18.42
QD-PMCD1002	1/8"	1/8" MT Pipe thread body w/Shut-off	15.68
QD-PMCD1004	1/8"	1/4" MT Pipe thread body w/Shut-off	15.68
QD-PLCD10006	1/4"	3/8" MT Pipe thread body w/Shut-off	12.91



MALE PIPE THREAD POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC100212	1/8"	1/8" MT Pipe thread body Straight-thru	15.52
QD-PMC100412	1/8"	1/4" MT Pipe thread body Straight-thru	15.76
QD-PLC1000412	1/4"	1/4" MT Pipe thread body Straight-thru	19.28
QD-PLC1000612	1/4"	3/8" MT Pipe thread body Straight-thru	17.66
QD-PMCD100212	1/8"	1/8" MT Pipe thread body w/Shut-off	17.33
QD-PMCD100412	1/8"	1/4" MT Pipe thread body w/Shut-off	17.93
QD-PLCD1000412	1/4"	3/8" MT Pipe thread body w/Shut-off	17.99
QD-PLCD1000612	1/4"	3/8" MT Pipe thread body w/Shut-off	19.19



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

LOW FLOW COUPLING INSERTS

IN-LINE HOSE BARB ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC1701	1/8"	1/16" Barb inline body Straight-thru	12.76
QD-PMC1702	1/8"	1/8" Barb inline body Straight-thru	8.66
QD-PMC1704	1/8"	1/4" Barb inline body Straight-thru	8.50
QD-PLC17004	1/4"	1/4" Barb inline body Straight-thru	8.99
QD-PLC17006	1/4"	3/8" Barb inline body Straight-thru	9.99
QD-PMCD1701	1/8"	1/16" Barb inline body w/Shut-off	18.76
QD-PMCD1702	1/8"	1/8" Barb inline body w/Shut-off	19.32
QD-PMCD1704	1/8"	1/4" Barb inline body w/Shut-off	19.32
QD-PLCD17004	1/4"	1/4" Barb inline body w/Shut-off	13.78
QD-PLCD17006	1/4"	3/8" Barb inline body w/Shut-off	13.66



IN-LINE HOSE BARB POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC170212	1/8"	1/16" Barb inline body Straight-thru	22.72
QD-PMC170412	1/8"	1/4" Barb inline body Straight-thru	23.54
QD-PLC1700412	1/4"	1/4" Barb inline body Straight-thru	17.92
QD-PLC1700612	1/4"	3/8" Barb inline body Straight-thru	18.77
QD-PMCD170112	1/8"	1/16" Barb inline body w/Shut-off	28.19
QD-PMCD170212	1/8"	1/8" Barb inline body w/Shut-off	22.72
QD-PMCD170412	1/8"	1/4" Barb inline body w/Shut-off	26.69
QD-PLCD1700612	1/4"	3/8" Barb inline body w/Shut-off	21.90



QUICK DISCONNECT FITTINGS

54

LOW FLOW COUPLING INSERTS

PANEL-MOUNT HOSE BARB ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC1601NA	1/8"	1/16" Barb panel mount body Straight-thru	17.49
QD-PMC1602NA	1/8"	1/8" Barb panel mount body Straight-thru	16.11
QD-PMC1604NA	1/8"	1/4" Barb panel mount body Straight-thru	16.16
QD-PLC16004NA	1/4"	1/4" Barb panel mount body Straight-thru	19.85
QD-PMCD1601NA	1/8"	1/16" Barb panel mount body w/Shut-off	25.34
QD-PMCD1602NA	1/8"	1/8" Barb panel mount body w/Shut-off	25.35
QD-PMCD1604NA	1/8"	1/4" Barb panel mount body w/Shut-off	28.34
QD-PLCD16004NA	1/4"	1/4" Barb panel mount body w/Shut-off	22.96
QD-PLCD16006NA	1/4"	3/8" Barb panel mount body w/Shut-off	22.96



PANEL-MOUNT HOSE BARB POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC160112NA	1/8"	1/16" Barb panel mount body Straight-thru	18.58
QD-PMC160212NA	1/8"	1/8" Barb panel mount body Straight-thru	18.73
QD-PMC16042NA	1/8"	1/4" Barb panel mount body Straight-thru	16.32
QD-PLC1600412NA	1/4"	1/4" Barb panel mount body Straight-thru	26.46
QD-PLC1600612NA	1/4"	3/8" Barb panel mount body Straight-thru	23.62
QD-PMCD160112NA	1/8"	1/16" Barb panel mount body w/Shut-off	20.63
QD-PMCD160212NA	1/8"	1/8" Barb panel mount body w/Shut-off	20.98
QD-PMCD160412NA	1/8"	1/4" Barb panel mount body w/Shut-off	21.52
QD-PLCD1600412NA	1/4"	1/4" Barb panel mount body w/Shut-off	26.46
QD-PLCD1600612NA	1/4"	3/8" Barb panel mount body w/Shut-off	26.96



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

QUICK DISCONNECT FITTINGS

55

LOW FLOW COUPLING INSERTS

IN-LINE FERRULESS COMPRESSION FITTING ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC1304NA	1/8"	1/4" OD Inline body Straight-thru	14.42
QD-PLC13006NA	1/4"	3/8" OD Inline body Straight-thru	18.92
QD-PMCD1304NA	1/8"	1/4" OD Inline body w/Shut-off	16.74
QD-PLCD13006NA	1/4"	3/8" OD Inline body w/Shut-off	19.99

IN-LINE FERRULESS COMPRESSION FITTING POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC130412NA	1/8"	1/4" OD Inline body Straight-thru	15.74
QD-PLC1300612NA	1/4"	3/8" OD Inline body Straight-thru	17.24
QD-PMCD130412NA	1/8"	1/4" OD Inline body w/Shut-off	28.19
QD-PLCD1300612NA	1/4"	3/8" OD Inline body w/Shut-off	30.91



PANEL-MOUNT FERRULESS COMPRESSION FITTING ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC1204NA	1/8"	1/4" OD Panel mount body Straight-thru	18.49
QD-PLC12006NA	1/4"	3/8" OD Panel mount body Straight-thru	29.47
QD-PMCD1204NA	1/8"	1/4" OD Panel mount body w/Shut-off	18.74
QD-PLCD12006NA	1/4"	3/8" OD Panel mount body w/Shut-off	26.41



PANEL-MOUNT FERRULESS COMPRESSION FITTING POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC120412NA	1/8"	1/4" OD Panel mount body Straight-thru	22.49
QD-PLC1200612NA	1/4"	3/8" OD Panel mount body Straight-thru	20.96
QD-PMCD120412NA	1/8"	1/4" OD Panel mount body w/Shut-off	23.98
QD-PLCD1200612NA	1/4"	3/8" OD Panel mount body w/Shut-off	32.43

MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

LOW FLOW COUPLING INSERTS

STRAIGHT HOSE BARB ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC2201	1/8"	1/16" Barb CPLG INS Straight-thru	8.38
QD-PMC2202	1/8"	1/8" Barb CPLG INS Straight-thru	5.38
QD-PMC2204	1/8"	1/4" Barb CPLG INS Straight-thru	5.81
QD-PMC22004	1/8"	1/4" Barb CPLG INS Straight-thru	5.66
QD-PLC22006	1/4"	3/8" Barb CPLG INS Straight-thru	3.06
QD-PMCD2201	1/8"	1/16" Barb CPLG INS w/Shut-off	10.28
QD-PMCD2202	1/8"	1/8" Barb CPLG INS w/Shut-off	9.26
QD-PMCD2204	1/8"	1/8" Barb CPLG INS w/Shut-off	9.26
QD-PLCD22004	1/8"	1/4" Barb CPLG INS w/Shut-off	9.48
QD-PLCD22006	1/8"	3/8" Barb CPLG INS w/Shut-off	9.48



STRAIGHT HOSE BARB POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC220112	1/8"	1/16" Barb valve hose Straight-thru	7.47
QD-PMC220212	1/8"	1/8" Barb valve hose Straight-thru	5.95
QD-PMC220412	1/8"	1/4" Barb valve hose Straight-thru	6.21
QD-PLC2200412	1/4"	1/4" Barb valve hose Straight-thru	6.21
QD-PLC2200612	1/4"	3/8" Barb valve hose Straight-thru	6.21
QD-PMCD220112	1/8"	1/16" Barb valve hose w/Shut-off	20.53
QD-PMCD220212	1/8"	1/8" Barb valve hose w/Shut-off	19.45
QD-PMCD220412	1/8"	1/4" Barb valve hose w/Shut-off	25.85
QD-PLCD2200412	1/4"	1/4" Barb valve hose w/Shut-off	15.98
QD-PLCD2200612	1/4"	3/8" Barb valve hose w/Shut-off	17.41

MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

LOW FLOW COUPLING INSERTS

STRAIGHT FERRULESS COMPRESSION FITTING ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC2004NA	1/8"	1/4" OD Compression Straight-thru	10.42
QD-PMC2006NA	1/8"	3/8" OD Compression Straight-thru	10.88
QD-PLC20006NA	1/4"	3/8" OD Compression Straight-thru	8.93
QD-PMCD2004NA	1/8"	1/4" OD Compression w/Shut-off	10.43
QD-PMCD2006NA	1/8"	3/8" OD Compression w/Shut-off	11.78
QD-PLCD20006NA	1/4"	3/8" OD Compression w/Shut-off	12.45



STRAIGHT FERRULESS COMPRESSION FITTING POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC200412NA	1/8"	1/4" OD Compression Straight-thru	10.20
QD-PLC2000612NA	1/4"	3/8" OD Compression Straight-thru	11.98
QD-PMCD200412NA	1/8"	1/4" OD Compression w/Shut-off	26.88
QD-PLCD2000612NA	1/4"	3/8" OD Compression w/Shut-off	22.19



PANEL MOUNT FERRULESS COMPRESSION FITTING ACETAL

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC4004NA	1/8"	1/4" OD Panel mount compression Straight-thru	11.12
QD-PLC40004NA	1/4"	1/4" OD Panel mount compression Straight-thru	15.46
QD-PLC40006NA	1/4"	3/8" OD Panel mount compression Straight-thru	16.21
QD-PMCD4004NA	1/8"	1/4" OD Panel mount compression w/Shut-off	15.73
QD-PLCD40004NA	1/4"	1/4" OD Panel mount compression w/Shut-off	17.76
QD-PLCD40006NA	1/4"	3/8" OD Panel mount compression w/Shut-off	18.63



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

LOW FLOW COUPLING INSERTS

ELBOW HOSE BARB ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC2302	1/8"	1/8" ID Barbed elbow Straight-thru	6.47
QD-PMC2304	1/8"	1/4" ID Barbed elbow Straight-thru	6.56
QD-PLC23004	1/4"	1/4" ID Barbed elbow Straight-thru	7.98
QD-PLC23006	1/4"	3/8" ID Barbed elbow Straight-thru	6.48
QD-PMCD2302	1/8"	1/8" ID Barbed elbow w/Shut-off	12.63
QD-PMCD2304	1/8"	1/4" ID Barbed elbow w/Shut-off	12.63
QD-PLCD23004	1/4"	1/4" ID Barbed elbow w/Shut-off	12.23
QD-PLCD23006	1/4"	3/8" ID Barbed elbow w/Shut-off	12.23



ELBOW HOSE BARB POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC230212	1/8"	1/8" ID Barbed elbow Straight-thru	8.57
QD-PMC230412	1/8"	1/4" ID Barbed elbow Straight-thru	8.57
QD-PLC2300412	1/4"	1/4" ID Barbed elbow Straight-thru	9.96
QD-PLC2300612	1/4"	3/8" ID Barbed elbow Straight-thru	9.96
QD-PMCD230212	1/8"	1/8" ID Barbed elbow w/Shut-off	15.81
QD-PMCD230412	1/8"	1/4" ID Barbed elbow w/Shut-off	15.81
QD-PLCD2300412	1/4"	1/4" ID Barbed elbow w/Shut-off	14.85
QD-PLCD2300612	1/4"	3/8" ID Barbed elbow w/Shut-off	15.52



PANEL MOUNT HOSE BARB ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC4202NA	1/8"	1/8" ID Barb panel-mount Straight-thru	11.83
QD-PMC4204NA	1/8"	1/4" ID Barb panel-mount Straight-thru	11.96
QD-PLC42004NA	1/4"	1/4" ID Barb panel-mount Straight-thru	14.93
QD-PLCD42004NA	1/4"	1/4" ID Barb panel-mount w/Shut-off	16.36
QD-PLC42006NA	1/4"	3/8" ID Barb panel-mount Straight-thru	11.67



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

QUICK DISCONNECT FITTINGS

59

LOW FLOW COUPLING INSERTS

ELBOW FERRULESS COMPRESSION FITTING ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC2104NA	1/8"	1/4" OD Elbow Straight-thru	9.69
QD-PMC2106NA	1/8"	3/8" OD Elbow Straight-thru	9.66
QD-PMCD2104NA	1/8"	1/4" OD Elbow w/Shut-off	17.67
QD-PLCD21006NA	1/4"	3/8" OD Elbow w/Shut-off	17.47



ELBOW FERRULESS COMPRESSION FITTING POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PLC2100612NA	1/4"	3/8" OD PP Elbow Straight-thru	14.23
QD-PMCD210412NA	1/8"	1/4" OD PP Elbow w/Shut-off	22.93
QD-PLCD2100612NA	1/4"	3/8" OD PP Elbow w/Shut-off	22.63

PIPE THREAD ACETAL (white)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC2402	1/8"	1/8" Male pipe thread Straight-thru	7.69
QD-PMC2404	1/8"	1/4" Male pipe thread Straight-thru	7.69
QD-PLC24004	1/4"	1/4" Male pipe thread Straight-thru	7.62
QD-PLC24006	1/4"	3/8" Male pipe thread Straight-thru	7.99
QD-PMCD2402	1/8"	1/8" Male pipe thread w/Shut-off	13.76
QD-PMCD2404	1/8"	1/4" Male pipe thread w/Shut-off	14.28



PIPE THREAD POLYPRO (beige)

PART NUMBER	FLOW	PRODUCT DESCRIPTION	PRICE \$
QD-PMC240212	1/8"	1/8" PP Male pipe thread Straight-thru	9.53
QD-PLC2400412	1/4"	1/8" PP Male pipe thread Straight-thru	9.35
QD-PMCD240212	1/8"	1/8" PP Male pipe thread w/Shut-off	15.93
QD-PLCD2400412	1/4"	1/4" PP Male pipe thread w/Shut-off	17.96



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

PARKER MALE CONNECTOR

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4MC2	1/4" Tube Compression X 1/8" Male NPT	3.84
P4MC4	1/4" Tube Compression X 1/4" Male NPT	3.99
P5MC2	5/16" Tube Compression X 1/8" Male NPT	8.94
P5MC4	5/16" Tube Compression X 1/4" Male NPT	4.48
P6MC2	3/8" Tube Compression X 1/8" Male NPT	5.22
P6MC4	3/8" Tube Compression X 1/4" Male NPT	4.84
P6MC6	3/8" Tube Compression X 3/8" Male NPT	5.48
P6MC8	3/8" Tube Compression X 1/2" Male NPT	5.54
P6MC12	3/8" Tube Compression X 1/8" Male NPT	8.97
P8MC2	1/2" Tube Compression X 1/8" Male NPT	8.99
P8MC4	1/2" Tube Compression X 1/4" Male NPT	5.99
P8MC6	1/2" Tube Compression X 3/8" Male NPT	6.46
P8MC8	1/2" Tube Compression X 1/2" Male NPT	6.90
P8MC12	1/2" Tube Compression X 3/4" Male NPT	7.97
P10MC2	5/8" Tube Compression X 1/8" Male NPT	9.67
P10MC4	5/8" Tube Compression X 1/4" Male NPT	9.66
P10MC6	5/8" Tube Compression X 3/8" Male NPT	8.83
P10MC8	5/8" Tube Compression X 1/2" Male NPT	8.64
P10MC12	5/8" Tube Compression X 3/4" Male NPT	9.88



PARKER FEMALE CONNECTOR

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4FC2	1/4" Tube Compression X 1/8" Female NPT	8.34
P4FC4	1/4" Tube Compression X 1/4" Female NPT	6.48
P6FC4	3/8" Tube Compression X 1/4" Female NPT	6.08
P6FC6	3/8" Tube Compression X 3/8" Female NPT	9.34
P6FC8	3/8" Tube Compression X 1/2" Female NPT	8.82
P8FC6	1/2" Tube Compression X 3/8" Female NPT	13.40
P8FC8	1/2" Tube Compression X 1/2" Female NPT	11.30
P10FC8	5/8" Tube Compression X 1/2" Female NPT	9.66



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

PARKER FITTINGS

61

PARKER MALE BRANCH TEE

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4MT2	1/4" Tube Compression X 1/8" Male NPT	10.80
P4MT4	1/4" Tube Compression X 1/4" Male NPT	9.54
P5MT2	5/16" Tube Compression X 1/8" Male NPT	14.98
P5MT4	5/16" Tube Compression X 1/4" Male NPT	16.32
P6MT4	3/8" Tube Compression X 1/4" Male NPT	10.68
P6MT6	3/8" Tube Compression X 3/8" Male NPT	12.74
P8MT6	1/2" Tube Compression X 3/8" Male NPT	18.54
P8MT8	1/2" Tube Compression X 1/2" Male NPT	19.66
P10MT8	5/8" Tube Compression X 1/2" Male NPT	14.64



PARKER MALE ELBOW

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4ME2	1/4" Tube Compression X 1/8" Male NPT	5.54
P4ME4	1/4" Tube Compression X 1/4" Male NPT	5.06
P4ME6	1/4" Tube Compression X 3/8" Male NPT	6.08
P5ME2	5/16" Tube Compression X 1/8" Male NPT	6.99
P5ME4	5/16" Tube Compression X 1/4" Male NPT	5.43
P5ME6	5/16" Tube Compression X 3/8" Male NPT	10.48
P6ME4	3/8" Tube Compression X 1/4" Male NPT	6.80
P6ME6	3/8" Tube Compression X 3/8" Male NPT	7.18
P6ME8	3/8" Tube Compression X 1/2" Male NPT	5.29
P8ME12	1/2" Tube Compression X 3/4" Male NPT	10.88
P8ME8	1/2" Tube Compression X 1/2" Male NPT	8.98
P10ME8	5/8" Tube Compression X 1/2" Male NPT	9.30



PARKER FEMALE ELBOW

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4FE2	1/4" Tube Compression X 1/8" Female NPT	7.58
P4FE4	1/4" Tube Compression X 1/4" Female NPT	7.98
P5FE2	1/4" Tube Compression X 5/16" Female NPT	13.76
P6FE4	3/8" Tube Compression X 1/4" Female NPT	8.88
P6FE6	3/8" Tube Compression X 3/8" Female NPT	6.92
P8FE6	1/2" Tube Compression X 3/8" Female NPT	9.68
P8FE8	1/2" Tube Compression X 1/2" Female NPT	10.80
P10FE8	5/8" Tube Compression X 5/8" Female NPT	9.70



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

PARKER TEE UNION

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4TU4	1/4" X 1/4" Tube Compression Tee	11.68
P5TU5	5/16" X 5/16" Tube Compression Tee	14.98
P6TU6	3/8" X 3/8" Tube Compression Tee	12.04
P6TU10	3/8" X 5/8" Tube Compression Tee	18.82
P8TU6	1/2" X 3/8" Tube Compression Tee	12.99
P8TU8	1/2" X 1/2" Tube Compression Tee	13.17
P10TU6	5/8" X 3/8" Tube Compression Tee	19.99
P10TU10	5/8" X 5/8" Tube Compression Tee	19.99



PARKER UNION CONNECTOR

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4UC4	1/4" X 1/4" Tube Compression Union	5.84
P5UC4	5/16" X 1/4" Tube Compression Union	7.74
P5UC5	5/16" X 5/16" Tube Compression Union	7.30
P6UC4	3/8" X 1/4" Tube Compression Union	7.44
P6UC5	3/8" X 5/16" Tube Compression Union	8.46
P6UC6	3/8" X 3/8" Tube Compression Union	7.26
P8UC6	1/2" X 3/8" Tube Compression Union	8.84
P8UC8	1/2" X 1/2" Tube Compression Union	8.96
P10UC6	5/8" X 3/8" Tube Compression Union	10.42
P10UC8	5/8" X 1/2" Tube Compression Union	10.86
P10UC10	5/8" X 5/8" Tube Compression Union	11.58



PARKER ELBOW UNION

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
P4EU4	1/4" X 1/4" Tube Comp Elbow Union	8.24
P5EU4	5/16" X 1/4" Tube Comp Elbow Union	14.14
P5EU5	5/16" X 5/16" Tube Comp Elbow Union	14.24
P6EU4	3/8" X 1/4" Tube Comp Elbow Union	12.60
P6EU5	3/8" X 5/16" Tube Comp Elbow Union	12.77
P6EU6	3/8" X 3/8" Tube Comp Elbow Union	10.54
P8EU6	1/2" X 3/8" Tube Comp Elbow Union	10.99
P8EU8	1/2" X 1/2" Tube Comp Elbow Union	12.86
P10EU10	5/8" X 5/8" Tube Comp Elbow Union	13.99



MUR-LOK COMPRESSION TEE

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-A0630405	3/8" Compression X 1/4" Female NPT	19.75
ML-A0630425	3/8" Compression X 1/4" ML	19.75
ML-A0630625	3/8" Compression X 3/8" ML	19.75
ML-A0630825	3/8" Compression X 1/2" ML	19.75
ML-A0830405	1/2" Compression X 1/4" Female NPT	19.75
ML-A0830425	1/2" Compression X 1/4" ML	19.75
ML-A0830625	1/2" Compression X 3/8" ML	19.75
ML-A0830825	1/2" Compression X 1/2" ML	19.75
ML-A1030405	5/8" Compression X 1/4" Female NPT	19.75
ML-A1030425	5/8" Compression X 1/4" ML	19.75
ML-A1030625	5/8" Compression X 3/8" ML	19.75
ML-A1030825	5/8" Compression X 1/2" ML	19.75



MUR-LOK IN-LINE BALL VALVES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-G0400405	1/4" Female NPT X 1/4" Female NPT	14.25
ML-G0400415	1/4" Female NPT X 1/4" Male NPT	14.25
ML-G0420405	1/4" ML X 1/4" Female NPT	14.25
ML-G0420415	1/4" ML X 1/4" Male NPT	14.25
ML-G0420425	1/4" ML X 1/4" ML	14.25
ML-G0420625	1/4" ML X 3/8" ML	14.25
ML-G0420825	1/4" ML X 1/2" ML	14.25
ML-G0600605	3/8" Female NPT X 3/8" Female NPT	14.25
ML-G0600615	3/8" Female NPT X 3/8" Male NPT	14.25
ML-G0620405	3/8" ML X 1/4" Female NPT	14.25
ML-G0620415	3/8" ML X 1/4" Male NPT	14.25
ML-G0620615	3/8" ML X 3/8" Male NPT	14.25
ML-G0620625	3/8" ML X 3/8" ML	14.25
ML-G0620825	3/8" ML X 1/2" ML	14.25
ML-G0820615	1/2" ML X 3/8" Male NPT	14.25
ML-G0820825	1/2" ML X 1/2" ML	14.25



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

MUR-LOK FITTINGS

64

MURLOK MALE NPTF CONNECTORS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-O0420216	1/4" ML X 1/8" Male NPT	1.93
ML-O0420416	1/4" ML X 1/4" Male NPT	1.96
ML-O0420616	1/4" ML X 3/8" Male NPT	2.96
ML-O0520416	5/16" ML X 1/4" Male NPT	1.98
ML-O0620216	3/8" ML X 1/8" Male NPT	3.96
ML-O0620416	3/8" ML X 1/4" Male NPT	2.91
ML-O0620616	3/8" ML X 3/8" Male NPT	3.54
ML-O0620816	3/8" ML X 1/2" Male NPT	5.53
ML-O0820616	1/2" ML X 3/8" Male NPT	4.36
ML-O0820816	1/2" ML X 1/2" Male NPT	4.88



MUR-LOK FIXED NPTF ELBOWS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-P0420206	1/4" ML X 1/8" Female NPT	3.45
ML-P0420216	1/4" ML X 1/8" Male NPT	2.81
ML-P0420416	1/4" ML X 1/4" Male NPT	2.98
ML-P0420616	1/4" ML X 3/8" Male NPT	3.46
ML-P0520416	5/16" ML X 1/4" Male NPT	3.63
ML-P0620416	3/8" ML X 1/4" Male NPT	3.93
ML-P0620616	3/8" ML X 3/8" Male NPT	4.48



MUR-LOK UNION & REDUCING ELBOWS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-Q0420426	1/4" ML X 1/4" ML	4.46
ML-Q0520426	5/16" ML X 1/4" ML	6.14
ML-Q0520526	5/16" ML X 5/16" ML	4.84
ML-Q0620426	3/8" ML X 1/4" ML	8.09
ML-Q0620526	3/8" ML X 5/16" ML	8.09
ML-Q0620626	3/8" ML X 3/8" ML	6.71
ML-Q0820626	1/2" ML X 3/8" ML	8.43
ML-Q0820826	1/2" ML X 1/2" ML	8.41



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

MUR-LOK FITTINGS

65

MUR-LOK REDUCERS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-S0420566	1/4" ML X 5/16" Stem	2.44
ML-S0420666	1/4" ML X 3/8" Stem	2.93
ML-S0520666	5/16" ML X 3/8" Stem	3.24
ML-S0620866	3/8" ML X 1/2" Stem	5.86



MUR-LOK NPTF STEM ADAPTERS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-S0460216	1/4" Stem X 1/8" Male NPT	1.98
ML-S0460416	1/4" Stem X 1/4" Male NPT	1.98
ML-S0660416	3/8" Stem X 1/4" Male NPT	2.48
ML-S0660616	3/8" Stem X 3/8" Male NPT	2.63
ML-S0860816	1/2" Stem X 1/2" Male NPT	2.85



MUR-LOK TUBE TO HOSE STEM ELBOWS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-S0460448	1/4" Stem X 1/4" BARB	2.96



MUR-LOK TUBE TO HOSE STEMS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-S0460446	1/4" Stem X 1/4" BARB	2.98
ML-S0660646	3/8" Stem X 3/8" BARB	3.18
ML-S0660846	3/8" Stem X 1/2" BARB	3.98
ML-S0860846	1/2" Stem X 1/2" BARB	3.80



MEDICAL SOLUTIONS INTERNATIONAL
5646 MERRIAM DRIVE
MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
913.438.9700
FAX 913.438.9701

MUR-LOK FITTINGS-PLUG IN ELBOWS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-T0420466	1/4" ML X 1/4" Stem	3.93
ML-T0420666	1/4" ML X 3/8" Stem	4.63
ML-T0520566	5/16" ML X 5/16" Stem	4.46
ML-T0620666	3/8" ML X 3/8" Stem	5.32
ML-T0820866	1/2" ML X 1/2" Stem	5.94



MUR-LOK FITTINGS-UNION TEES

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-V0420426	1/4" ML X 1/4" ML X 1/4" ML	4.68
ML-V0420526	1/4" ML X 1/4" ML X 5/16" ML	6.18
ML-V0520526	5/16" ML X 5/16" ML X 5/16" ML	4.82
ML-V0620426	3/8" ML X 3/8" ML X 1/4" ML	6.37
ML-V0620626	3/8" ML X 3/8" ML X 3/8" ML	5.79
ML-V0820426	1/2" ML X 1/2" ML X 1/4" ML	10.96
ML-V0820626	1/2" ML X 1/2" ML X 3/8" ML	7.24
ML-V0820826	1/2" ML X 1/2" ML X 1/2" ML	7.56
ML-V0820866	1/2" ML X 1/2" ML X 1/2" Stem	6.56



MUR-LOK FITTINGS-FEMALE CONNECTORS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
ML-U0420206	1/4" ML X 1/8" Female NPTF	4.72
ML-U0420406	1/4" ML X 1/4" Female NPTF	6.56
ML-U0620406	3/8" ML X 1/4" Female NPTF	5.96



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

NYLON GARDEN HOSE FITTINGS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
GH-TVW7	"O" RING for Garden Hose Fitting	0.50
GH-TA772	3/4" Male Garden Hose X 1/4" Male NPT	3.45
GH-TA774	3/4" Male Garden Hose X 1/2" Male NPT	3.95
GH-TA776	3/4" Male Garden Hose X 3/4" Male NPT	4.05
GH-TA792	3/4" Male Garden Hose X 1/4" Female NPT	5.10
GH-TA793	3/4" Female Garden Hose X 3/8" Male NPT	5.10
GH-TCF72	3/4" Female Garden Hose X 1/4" Barb	7.61
GH-TCF73	3/4" Female Garden Hose X 3/8" Barb	7.61
GH-TCF76	3/4" Female Garden Hose X 3/4" Barb	8.79
GH-TCM72	3/4" Male Garden Hose X 1/4" Barb	4.93
GH-TCM73	3/4" Male Garden Hose X 3/8" Barb	3.95
GH-TTA-776	3/4" Female Garden Hose X 3/4" Male NPT	6.72



NYLON WALL STATION DRAIN BARBS

PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
W-ATE2	Wallstation Drain Barb 1" Male NPT X 3/4" Barb Elbow	3.70
W-ATA1	Wallstation Drain Barb 1" Male NPT X 3/4" Barb Straight	3.10



MEDICAL SOLUTIONS INTERNATIONAL
 5646 MERRIAM DRIVE
 MERRIAM, KS 66203
www.medicalsolutionskc.com

800.326.5275
 913.438.9700
 FAX 913.438.9701

TANKS		
PART NUMBER	PRODUCT DESCRIPTION	PRICE \$
T-SSTA-120	120 Gallon stackable tank	QUOTE
T-ST220	220 Gallon stackable tank w/ colored base	QUOTE
T-ST225	225 Gallon stackable tank w/colored base	QUOTE
T-SSTA-240	240 Gallon stackable tank	QUOTE
T-ST330	330 Gallon stackable tank w/colored base	QUOTE
T-10250-B	250 Gallon Dish Bottom w/ Blue Base	QUOTE
T-10250-O	250 Gallon Dish Bottom w/ Orange Base	QUOTE
T-10250-R	250 Gallon Dish Bottom w/ Red Base	QUOTE
T-10250-Y	250 Gallon Dish Bottom w/ Yellow Base	QUOTE
T-10300	300 Gallon storage tank	QUOTE
T-10500-FB	500 Gallon storage tank flat bottom	QUOTE
T-10500-CB	500 Gallon storage tank conical bottom	QUOTE
T-VC-200	Ventilation Cap for square stackable tank	24.99



ALL FORMS OF TANKS ARE AVAILABLE, WE ARE PLEASED TO QUOTE CUSTOM ORDERS.



ORDER FORM

DATE _____

PO# _____

BILL TO:

SHIP TO:

QUANTITY	ITEM#	DESCRIPTION	COST
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

MEDICAL SOLUTIONS INTERNATIONAL 5646 MERRIAM DRIVE, MERRIAM, KANSAS 66203
800-326-5275, 913-438-9700, FAX 913-438-9701

sales@medicalsolutionskc.com



Medical Solutions

I N T E R N A T I O N A L

5646 Merriam Drive Merriam, Kansas 66203
913-438-9700 Fax 913-438-9701 800-326-5275
•www.medicalsolutionskc.com